



Kiran Surgical is a technically driven partnership firm active in the commerce of Medical Disposable Products. Established in the year 2007, we, Kiran Surgical are an illustrious manufacturer and exporter of Medical Products. Our products are widely applied in varied medical organizations, hospitals, clinics and nursing homes. Due to the superior quality of our products and large product line, we have gained strong credibility in the market. Our firm has a contemporary infrastructure facility that helps us in exporting our products in bulk. This facility is equipped with the essential machines and tools required for the supplying and exporting process. Our range is offered using premium grade raw material, which is sourced from the trustworthy vendors of the market. Stringent quality testing measures are adopted to ensure that the offered assortment is flawless. Due to our immaculate range, we have achieved tremendous client satisfaction from different nook and corners of the world.

www.kiransurgicalsystems.com
www.kiranmedicareequipments.com
www.kisusy.com, www.d4h.co.in, www.d4s.co.in



Delhi : A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092

Gujarat : 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda – 390 010,

Andhra Pradesh : No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad – 500 034

Tamil Nadu : No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore - 641 021

Kerala : 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621

Helpline : +91 9751725000, +91 9751625000, + 91 9751425000, + 91 9751215000, + 91 9751245000

Email : kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com

Website : www.kirangroupofcompanies.com, www.kiran.biz, www.kiranshop.in, www.kiranindia.in, www.kisusy.com, www.kiransurgicalsystems.com, www.kiranmedicareequipments.com, www.kiranpump.com, www.d4h.co.in, www.kiransolarsystems.com, www.kiranled.com, www.kiranincinerationsystems.com, www.d4s.co.in, www.kiranequipments.com, www.kirancompressors.com, www.kiranappliances.com, www.kiranairventilators.com, www.kiranfrpproducts.com, www.kiranpebsystems.com, www.kiranpanels.com, www.kiranfabssystems.com

Products

Kiran Latex Surgical Gloves

Product Overview

Manufactured with Non-Toxic Chemical and Special quality Latex In Clean Hygienic and dust free environment. Pre-Powdered with Bio-Absorbable Cornstarch to prevent granulosa with extensive leaching in hot flowing water (80 ft.) for reducing non-latex residues to prevent allergic reactions. Sterilized by E.T.O. Uncompromised inspection and attention to each glove for added safety for Surgeon and Patients. High tensile Strength - Meets / Exceeds international Standards. Anatomically shaped for comfort, textured surface for excellent instrument grip.

- Size: 6", 6.5", 7", 7.5" and 8" Pair. (Individually packed Latex Gloves).
- Carton Packing: 50 pairs in a inner box, 1000 pairs in a carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)
- Export Packing: 180 Cartons in a 20' Container, (1,80,000 Pairs / 20' Container)

Kiran Latex Examination Gloves

Product Overview

Available with us, is an assortment of latex examination gloves which is made of quality material. These gloves are especially designed for efficient performance in various surgical and medical applications. Owing to comfort, durability and skin friendliness, these are in high demand. Our products are available in variegated sizes so as to meet the different requirements of the clients.

- Size: Small, Medium, Large and Ex-Large
- Carton Packing: 100 pieces in a inner box, 2500 Pieces in a carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)
- Export Packing: 180 Cartons in a 20' Container, (4,50,000 Pieces / 20' Container)

Plastic Examination Gloves

Product Overview

Protection of hand against contamination. Eg. While nursing or examining patient. Recommended for the use of gynecologist during per vagina or per rectum examination. Recommended for use in

Industries where safety and hygiene are of prime importance. Specifications :

Made from high quality plastic.

Safe & Hygienic

Fits either hand

Good barrier properties

- Carton Packing: 100 Pieces in a Bag, 10,000 Pieces in a Carton
- Carton Size: 30" (L) x 15" (W) x 10" (H)
- Export Packing: 360 Cartons in a 20' Container, (36,00,000 Pieces/20' Container)

Kiran HIV Kit

Product Overview

We are counted as one of the prime HIV Kit Manufacturer and Suppliers based in India. The HIV Kit offered by us consists of 7 items packed in Medicated Paper Pouch. The contents in HIV Kit are as under:

- | | | |
|-----------------------------|-------------|-----------|
| (1) Surgeon's Gown | - 40 gsm | - 1 Piece |
| (2) Face Mask | - Non-Woven | - 1 Piece |
| (3) Surgeon's Cap | | - 1 Piece |
| (4) Latex Surgical Gloves | | - 1 Pair |
| (5) Shoe Cover | - Non-Woven | - 1 Pair |
| (6) Goggles with Adjustment | | - 1 Piece |
| (7) Waste Carry Bag | | - 1 Piece |

- Carton Packing: Kit of 7 items in a Medicated Paper Pouch, 50 Kits in a carton.
- Carton Size: 26" (L) x 18" (W) x 21" (H)





Kiran Disposable Delivery Kit with Linen

No.	Item	Quantity / Kit
1	Surgical Blade	1 No.
2	Face Mask	3 Nos.
3	Cap	3 Nos.
4	Drape Sheet	2 Nos.
5	Sanitary Pad	5 Nos.
6	Perineal Support Pad	1 No.
7	Apron	3 Nos.
8	Latex Gloves	2 Nos.
9	Shoe Cover	3 Pairs
10	Umbilical Cord Clamp	3 Nos.
11	Baby Blanket	1 No.
12	Baby wrap towel	1 No.
13	Perineal Sheet	1 No.
14	Mucus Extractor	1 No.
15	Cotton Swab	2 Nos.

E.T.O. Sterile & Pyrogen Free. Also available non-sterile (Packed Clinically).

- Carton Packing: 15 Items in a Inner Box, 12 Boxes in a outer Carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)
- Export Packing: 180 Cartons in a 20' Container, (2,160 Pieces / 20' Container)

Kiran Oxygen Mask

Product Overview

Oxygen Mask is designed for convenient oxygen therapy; Mask is moulded from non-toxic medical grade PVC. It is provided with adjustable strap and nasal oxygen mask integrated Nose Bridge for proper positioning of mask on the nasal area. Other features includes:

- Swivel connector for convenience of attachment to the oxygen tube.
- 2 meter long multi channel tube to ensure continuous flow of oxygen.
- Proximal end of tube is fitted with soft funnel shape connector for easy connection to oxygen source.
- E.T.O. Sterile & Pyrogen Free. Also available non-sterile (Packed Clinically).

Size: Adult and Pediatric Colour: Green/Clear

- Carton Packing: Individually packed in a bag, 75 pieces in a carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)
- Export Packing: 180 Cartons in a 20' Container, (13,500 Pieces / 20' Container)

Kiran Urine Bag

Product Overview

Manufactured from Flexible and Ultra Soft Pvc Film Provides Non Return Valve to Prevent Back Flow of Urine Capacity - 2000 ML Tapered Connector Clear Transparent, Kink Resistant Pvc Tubing.

- Carton Packing: 10 pieces in a inner bag, 400 pieces in a carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)

Kiran Shoe Cover – Plastic

Product Overview

We offer a wide collection of PP shoes cover which is manufactured from high quality raw material. PP shoe covers are made of high grade plastic which makes it hygienic and highly durable. Our range of PP shoe covers are highly appreciated by our clients as it has got universal fitting and dust resistive. The main purpose of PP shoe cover is to provide safety against germs and keeps the work place dust free. PP shoe covers are available in blue color.

Plastic material used :

- HMHD Plastic
- LLDPE plastic

- Carton Packing : 50 Pairs in a plain PP bag, 3000 pairs in a carton
- Carton Size: 30" (L) x 16" (W) x 20" (H)



Products

Kiran Umbilical Cord Clamp

Product Overview

Offered by us is an Umbilical Cord Clamp, which is demanded in various medical processes. The range of our services is rendered using quality metal and material by the diligent professionals who are designed for using immediately after the birth. Besides these products are provided with double purpose security lock and protection against accidental re-opening. These products have special attributes for which these are highly demanded for the following attributes.

- Designed for clamping the umbilical cord immediately after birth
- Finger grip ensures safe & convenient handling.
- Provided with grooves all along the length.
- Manufactured from non toxic medical grade polymer.
- E.T.O. sterile & pyrogen free.
- Packing: Individually packed in peelable paper pouch.,
Inner Box of 200 Pieces
- Carton Packing: 200 pieces in a Inner Box, 2400 pieces in a carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)



Kiran Male External Catheter

Product Overview

Penile sheath/ external catheter, provided with self adhesive coated strip for proper fixation on to the penis. Proximal end is designed for easy connection to urine bag/leg bag, making it simple to use. Individually packed and ready for use.

Size: Small	(20 mm Diameter)
Medium	(25 mm Diameter)
Large	(30 mm Diameter)
Ex-Large	(40 mm Diameter)

- Packing: Individually packed in peelable paper pouch.,
Box of 50 Pieces
- Carton Packing: 50 pieces in a inner Box, 600 pieces in a carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)



Kiran Surgeon Gown

Product Overview

The Medical Surgeon Gown that we cater is made from high quality spun bond polymer non-woven fabric and spun lace fabrics making them light in weight. Surgeon Gown is designed to prevent the body and clothing from contact with dust and liquid. Besides, we are one of the distinguished Disposable Surgeon Gown Exporters and Suppliers in India.

Features of Disposable Non Woven Gown

Universal fittings, prevents body and clothing from contact with dust and liquid. Most protective and economic, stops migration of skin cell and virtually lint free. Soft light weight, comfortable and provides opacity, stronger fabric resists tearing.

Technical Specification of Disposable Non-Woven Gown

Knitted Cuff

5 Thread stitching

Velcro for locking

40 GSM

120 CM. Length

- Carton Packing : Individually packed in plain PP bag, 150 pieces in a carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)
- Export Packing: 180 Cartons in a 20' Container, (27,000 Pieces / 20' Container)





Kiran Guedel Airway (Oralpharyngeal Airway)

Product Overview

Owing to our experience, we have been successful in catering to the requirements of our esteem clients by offering quality Guedel Airway. These products are made of color coded connector which prevents occlusion and tough biting. Apart from that, our products are pyrogen free and E.T.O sterilized as well. Further these products have special attributes for which they are highly demanded in the market:

- Manufactured from EVA to achieve extra soft surface finish
- The integral hard bite block avoids airway occlusion and biting of the tongue
- Colour coded connector for easy size identification.
- Individually packed in a Medicated Paper pouch.

Number	Use for	Size	Colour
0	New Born	40 mm	Pink
1	Infant	50 mm	Blue
2	Small Child	60 mm	Black
3	Child	70 mm	White
4	Small Adult	80 mm	Green
5	Medium Adult	90 mm	Yellow
6	Large Adult	100 mm	Red

- Carton Packing: Individually Packed in Medicated Paper Pouch, Size: 00 to 3 (50 pieces in a Inner Box), 4 and 5 (40 Pieces in a Inner Box)
- 12 Inner boxes in a master carton.
- Carton Size: 26" (L) x 18" (W) x 21" (H)

Kiran Infant Mucus Extractor

Product Overview

Our company caters to a wide market for Infant Mucus Extractor, which are specially designed to free respiration problems of infant. These products are transparent and the surfaces of these products have low friction. Apart from this, we are well known as one of the leading exporter of Disposable Infant Mucus Extractor. Following are the attributes of these products for which they are widely known : Specially designed for aspiration of secretion from oropharynx in newly born babies to ensure free respiration. Low friction surface catheter is provided with open end silk smooth round tip for trauma free insertion of the aspiration. Also suitable for obtaining mucus specimen for microbiological examination. Spare plug cap is provided to seal the container for safe transportation of specimen to the laboratory or aseptic disposal of container. • Packing: Individually packed in peelable paper pouch., Inner Box of 40 Pieces

- Carton Packing: 40 pieces in a inner Box, 480 pieces in a carton
- Carton Size: 26" (L) x 18" (W) x 21" (H)

Kiran Foleys Balloon Catheter – 2 way

Product Overview

Kiran Surgical is a well acclaimed supplier of high quality Foley Balloon Catheter. Foley Balloon Catheter is a Urinary tract catheter used for short or long term urine drainage. The Foley Balloon Catheter, we offer is made from natural and premium latex rubber.

Features

Siliconized smooth surface for non traumatic intubation

Distal end is coned with super smooth eyes proportional to each other to provide non-traumatic intubation. Provided with ultra thin, highly elastic Balloon. Hard Non Return valve ensures trouble free inflation and deflation of balloon.

Size: FG: 10, 12, 14, 16, 18, 20, 22 and 24.

- Carton Packing: Sterile Individually Packed in peelable Blister pack, 10 Pieces in a Inner Box, 600 Pieces / Carton.
- Carton Size: 26" (L) x 12.5" (W) x 21" (H)



Kiran Face Mask

Product Overview

We are one of the renowned Medical Surgical Face Mask Exporters based in the country. Surgical Face Mask is commonly used to protect the patient against airborne infective agents that are originated from the respiratory tract of the healthcare staff.

Features of Disposable Face Mask

Masks are primarily intended to protect the patient against airborne infective agents originating from the healthcare staff's respiratory tract. Kiran Surgical masks are the precise balance of bacteria filtration efficiency (BFE), Prevents skin irritation, pleated for full face coverage, Ultra sonic sealing & proper placement of folds. Cool & breathable, lint Free.

Technical Specification of Disposable Face Mask

All the above disposable face mask are made from spun bond polymer non-woven fabric, spun lace fabrics as the main raw material. This material is preferred for its light weight. Its gravity is lighter than hydrophobic PP fibres so it has excellent water repellence and air permeability resulting in providing comfort & dryness resulting in anti bacterial and anti fungal resistance.

Thickness: 2 Ply and 3 Ply

Type: Tie (Lace) and Elastic (Ear Loop)

Colour: Green, Blue and White

Packing: Individually Packed Face Mask

• Carton Packing : 50 Pieces in a inner box, 2200 Pieces in a carton

• Carton Size: 26" (L) x 18" (W) x 21" (H)

• Export Packing: 180 Cartons in a 20' Container, (3,96,000 Pieces / 20' Container)

Kiran Surgeon's Cap

Product Overview

We are offering Non Woven Surgeon Cap. Our medical surgeon cap is made from spun bond polymer non-woven fabric and spun lace fabrics. This cap prevents less movement of hair which contains a large number of micro organisms that can be transferred at the time of surgery. The cap is completely ultrasonically sealed

Features of disposable non-woven cap

The hair and scalp area frequently contain large numbers of micro organisms. If hair is kept still from movement then less bacteria is released and the most effective way to keep the hair still is to cover it entirely with surgical headwear. Kiran Surgical cap is completely ultrasonically sealed, giving a very high quality and hygienic finish. More economical and far superior than previous stitched caps.

Type: Tie (Lace) and Elastic

Colour: Green, Blue and White

Packing: Individually Packed Surgeon's Cap

• Carton Packing: 50 Pieces in a inner box, 2200 Pieces in a carton

• Carton Size: 26" (L) x 18" (W) x 21" (H)

Kiran Bouffant (Nurse) Cap

Product Overview

We are offering Non Woven Bouffant Cap. Our medical Bouffant Cap is made from spun bond polymer non-woven fabric and spun lace fabrics. This cap prevents less movement of hair which contains a large number of micro organisms that can be transferred at the time of surgery. The cap is completely ultrasonically sealed.

Features of disposable non-woven cap

The hair and scalp area frequently contain large numbers of micro organisms. If hair is kept still from movement then less bacteria is released and the most effective way to keep the hair still is to cover it entirely with surgical headwear. Kiran Surgical cap is completely ultrasonically sealed, giving a very high quality and hygienic finish. More economical and far superior than previous stitched caps.

Type: Round Elastic

Colour: Green, Blue and White

Packing: Loose 50 Pieces in a Plain PP bag

• Carton Packing: 50 Pieces in a inner bag, 5000 Pieces in a carton

• Carton Size: 27" (L) x 16" (W) x 19" (H)

Kiran Shoe Cover – Non Woven

Product Overview

We are offering prime quality non woven shoe cover which is made from durable and strong non woven fabric. Non woven shoe covers are mainly used in hospitals, nursing home, clinics, kitchens and other health sectors. These shoe cover are reliable form of protecting room areas, such as Operation Theater, etc. from dust and germs.

• Carton Packing : 25 Pairs in a plain PP bag, 1250 pairs in a carton

• Carton Size: 27" (L) x 18" (W) x 16" (H)

