



Kiran TECHNO SERVICES
PRIVATE LIMITED



KIRAN MEDICAL DISPOSABLES

www.kiransurgicalsystems.com
www.kiranmedicareequipments.com



DISPOSABLE SURGICAL PRODUCTS

We are the leading Manufacture, Exporter and Supplier of Disposable Surgical Product in India like Surgical Latex Gloves, Examination Gloves, Unsterile Gloves, Infusion Set, Blood Administration, Micro-Infusion Set, Scalp Vein Set, Urine Bag, Nelaton Catheter, Suction Catheter, Ryles Tube, Infant Feeding Tube, Pedia Drip Set, Face Mask, Disposable Cap, Shoe Covers, Disposable Apron, Gowns, Dipper, Condom Catheter, Foley Catheter, Karman Cannula, Oxygen Mask, Cord Clamp, Mucus Extractor etc.

INFUSION SET (MRS)



- Conforms to International Standards
- Special beveled extra strong spike
- Self sealing latex tube for extra medication & easy flushing
- Highly kink resistant tubing. Roller clamp for flow controller
- Double bag packing ensures sterility till product is actually used Needle 21G x 1 1/2" for smooth flow
- Disposable non toxic & pyrogen free, Drop size of 20 drops/ml.,
- Packing : 25 Pcs/Pkt.

MEASURED VOLUME SET (PEDIA DRIP)



- Burette type measured vol. of 110ml, with 10ml, overflow limit
- Efficient Roller controller for better flow control
- Specially designed floating auto shut off valve acts as floating indicator and automatically shut off the drain path when the chamber gets empty to prevent air trapping in the fluid line
- Separate plug for extra medication and continues changeover
- Drop size of 60 drops/ml, Cap.:110ml./150ml
- Disposable, non toxic & pyrogen free
- Pack.:10 pcs. Box

MEASURED VOLUME SET (PEDIA DRIP)



- Burette type measured vol. of 110ml, with 10ml, overflow limit
- Efficient Roller controller for better flow control
- Specially designed floating auto shut off valve acts as floating indicator and automatically shut off the drain path when the chamber gets empty to prevent air trapping in the fluid line
- Separate plug for extra medication and continues changeover
- Drop size of 60 drops/ml, Cap.:110ml./150ml
- Disposable, non toxic & pyrogen free
- Pack.:10 pcs. Box

FOLEY'S CATHETER (FOLEY CATH)



- Symmetrical large capacity balloon ensuring a straight tip for proper flow and good sphincter action thus preventing bladder leakage
- The eyes are burr free and bonded with silicone elastomer to eliminate the risk of injury and encrustation
- Valve ensures trouble free inflation and deflation of balloon
- Balloon capacity: 30ml
- SIZE: Fg: 8,10,12, 14, 16, 18, 20, 22 & 24

OXYGEN MASK (MEDI OXY MASK)



- Moulded face mask has adjustable elastic strap and integrated nose clip for proper positioning of mask on the mouth and nasal area
- Specially designed for convenient oxygen therapy
- Provided with 2 meter long star lumen tube to ensure continuous flow of Oxygen.
- The tube has funnel shaped connector at the proximal end to facilitate simple & safe connection with oxygen tube
- Sizes: Adult, Paediatric
- Packing: 10 Pcs/ Box

URINE COLLECTING BAG (URO-10)



- Manufactured from clinical grade PVC Sheet, Capacity 2000ml. 100cm long super-smooth kink resistant tubing is easy to milk and ensures unrestricted flow
- Specially designed hook facilities for carrying & handling of bag and holds the tube in upright position
- Efficient non return valve prevents back flow of bacteria migration
- Drainage outlet system Modified, enlarged and tapered connector facilitates aseptic catheter-urine bag connection

BLOOD TRANSFUSION (MRS B.T)



- Conforms to International Standard
- Transparent & flexible double drip chamber facilitates visual access and rapid adjustment of fluid level
- Extra strong bevel shaped spike for easy penetration. Highly kink resistant tubing
- Efficient roller control for accurate adjustment of infusion rate
- Self sealing latex bulb for extra medication & easy flushing
- Double packing ensures sterility till the product is actually used
- Needle 18 x 1 1/2" for smooth flow
- Drop size 20 drops/ml.,
- Disposable non toxic & pyrogen free
- Packing : 20 Pcs/Pkt.

INFANT MUCUS EXTRACTOR



- Specially designed for aspiration of secretion from Oropharynx in newly born babies to ensure free aspiration
- Clear transparent container permits immediate visual examination of the aspirate
- Also suitable for obtaining mucus specimen for micro biological examination
- Packing : 10 Pcs/ Box

DISPOSABLE CAP (MEDI CAP)



- Made from soft non-woven cloth
- Specially designed to protect the patient against dropping of hair/dandruff of the doctors and other attending staff
- Provided with stitched interlocking adaptable elastic for better grip on the forehead and comfortable fitting
- Individually packed
- Disposable, Sterile ready for use
- Packing: Box of 100 pcs

SUCTION CATHETER (SUCTION CATH)



- Manufactured from non toxic non irritant medical grade PVC
- Suction catheters are suitable for removal of secretion from mouth oropharynx trachea and bronchial tubes
- Distal eye is open with one lateraleye
- Frozen surface tubing for smooth intubation
- Colour Coded for instant size identification
- Length: 53 cms,
- Size: 8,10,12,14,16,18,20
- Packing: 10 Pcs/Pkt.

DISPOSABLE FACE MASK (MEDI MASK)



- Made from soft non Woven cloth
- Double layer mask offering excellent bacterial filtration efficiency
- Super smooth soft ultrasonically welded mask is most comfortable during use, Disposable, Sterile ready for use
- Individually packed.
- Packing: Box of 100 pcs

INFANT FEEDING TUBE



- Manufactured from non toxic PVC non irritant to delicate mucosa suitable for neonates and paediatric nutritional feeding
- Distal end is coned with two lateraleyes
- Proximal end is fitted with female luer mount for easy connection feeding funnel or syringe
- Low friction tube & super smooth tip ensures trauma free intubation
- Radio-opaque line is provided throughout its length for X-ray visualization. Length: 52 cm
- SIZE: FG 5, 6, 7, 8, 9, 10

LATEX EXAMINATION GLOVES



- Made from pure natural Latex for greater finger dexterity
- Non finger fatigue. Keep the hand cool & comfortable
- Laminated inner box of 100 pcs
- SIZE: Small, Medium & Large

LATEX SURGICAL GLOVES (Sterile)



- Made from pure white natural latex for greater finger dexterity
- Non finger fatigue. Keeps the hands cool and comfortable
- For single use only.
- Individually packed
- Disposable, Sterile ready for use, E.T.O. Sterile.
- Micro rough finish for better grip
- Curved fingers, anatomy shaped for extra comfort
- Sterile pre-powdered
- Size: 6, 6.5, 7, 7.5, 8.
- Packing: 25 Pcs/Box

NELATON CATHETER (NEL CATH)



- Manufactured from non-toxic non irritant PVC Compatible with catheter lubricants
- Nelaton Catheters are used for short term bladder catheterisation through urethra
- Distal end is coned for non traumatic introduction with two lateral eyes for efficient drainage
- Proximal end is fitted with universal funnel shaped connector for extension
- Frozen surface tubing for super smooth intubation
- Funnel shape connectors are colour coded for instant identification of sizes
- Length: 40 cms.
- Size: Fg: 8,10,12, 14, 16, 18

SCALP VEIN SET (MEDI VEIN SET)



- SV set is designed to provide rapid venous access with greater patient comfort during infusion Short beveled siliconised needle facilitates atraumatic cannulation
- Thin wall needle provides better flow rate, as in the same cannula more circumference is available for better flow of liquid
- Butterfly shaped wings facilitates easy handling & attachment with the skin
- Proximal end is provided with flexible female luer fitting
- Butterfly connected to soft non-toxic non irritant tube which does not kink or coil set
- Butterflies are colour coded for identification of needle size
- Size: G-18,19,20,21,22,23,24

MALE CATHETER (MALE CATH)



- Penile Sheath/External Catheter
- Manufactured from pure latex for soft and gentle feel
- Male catheter is specially designed for urine incontinence for day and night use in male patient
- Provided with self adhesive coated strip for proper fixing on to the penis
- Proximal end is designed for easy connection to urine bag, making it simple to use
- Size: Small, Medium & Large
- Packing: 25 Pcs/Box



www.kiranmedicareequipments.com

www.kiransurgicalsystems.com



Delhi : A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092

Gujarat : 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda – 390 010,

Andhra Pradesh : No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad – 500 034

Tamil Nadu : No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore - 641 021

Kerala : 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621

Helpline : +91 9751725000, +91 9751625000, +91 9751425000, +91 9751215000, +91 9751245000

Email : kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com

Website : www.kirangroupofcompanies.com, www.kiran.biz, www.kiranshop.in, www.kiranindia.in, www.kiransurgicalsystems.com, www.kiranmedicareequipments.com, www.kiranpump.com, www.kiransolarsystems.com, www.kiranled.com, www.kiranincinerationsystems.com, www.kiranequipments.com, www.kirancompressors.com, www.kiranappliances.com, www.kiranairventilators.com, www.kiranfrpproducts.com, www.kiranpebsystems.com, www.kiranpanels.com, www.kiranfabssystems.com