

KIRAN MEDICAL EQUIPMENTS



www.kiransurgicalsystems.com

Catheterization Lab

- 3 different software configuration available: DFG, MN12, MN25 all with DSA/angio cardio measurements
- E-MOTION: all C-arm movements can be motorized
- Dynamic anticollision system on all C-arm movements.
- DUAL COOLING SYSTEM: for immediate and effective heat removal. Innovative and exclusive liquid cooling system with liquid-to-air heat exchanger. At the nominal X-ray power the temperature of the X-ray assembly remains steady for an indefinite period of time.
- DUAL POWER SYSTEM: power reserve system allows functioning from any standard wall outlet. Power is stored and released according to necessity.
- DAP meter
- Full Dicom 3, USB, CD/DVD Interfacing
- Remote Control



Portable Ultrasound Systems

Being client centric organization, we are engaged in offering wide assortment of **Portable Ultrasound Systems**. We offer Portable Ultrasound Systems to our patrons as for diagnosis the patient's body internally. Our Portable Ultrasound Systems are used in labs and many more places and known for their accurate results. Moreover, we offer Portable Ultrasound Systems at affordable rates.

Features:

- | | |
|---|---|
| ◆ Elegant Image Quality for complete clinical diagnosis | ◆ Ideal Portable work station |
| ◆ 10" non - interlaced display | ◆ Dual Transducer connector - built in standard |
| ◆ Various types of high performance multi frequency transducers | |

Specifications:

- | | | | |
|-----------------------------|------------------------------|------------------------------|-----------------|
| ◆ 3.5 MHz convex probe | ◆ 5.0 MHz micro-convex probe | ◆ 6.5 MHz transvaginal probe | ◆ Video printer |
| ◆ 7.5 MHz transrectal probe | ◆ Application software (CD) | ◆ High speed USB cable | ◆ Battery |
| ◆ Power adapter | ◆ 7.5 MHz linear probe | | |



ABG Analyser OPTI CCA TS

- ◆ Innovative optical fluorescence and reflectance technology
- ◆ Only STAT portable device that measures THB and So2
- ◆ Provides fast accurate results for blood gases, acid base balance and electrolytes
- ◆ Lower cost per test as no external reagents are required
- ◆ Run samples of whole blood, plasma and serum on single use disposable cassettes
- ◆ Waste is safely sealed inside the cassette for clean, easy and safe biohazard disposal
- ◆ Room temperature storage for immediate use without warm-up
- ◆ No electrodes required provides maintenance-free operation
- ◆ Automatic sample aspiration for accurate results
- ◆ Results in less than 120 seconds
- ◆ Rechargeable battery backup of up to 8 hours
- ◆ Intuitive touch screen user interface
- ◆ Upgradable machine
- ◆ CE, UL, IVDD and FDA approved



MRI Machines

We are renowned amongst the customers for offering wide assortment of **MRI Machines**. Our products are used for diagnosis patient in labs, hospitals and many more places. Customers like our products for their long functional life and effectiveness. Moreover, we offer MRI Machines to our patrons in recommended specifications.


◆ Tesla	0.3 Tesla	◆ Magnet	Permanent	◆ Magnet Type	Open Type
◆ Gradient Strength	21mt/m	◆ Slew Rate	55mt/m/Sec	◆ Channel	4 Ch
◆ Gap	42.5 cm V	◆ Monitor type	CRT (LCD)	◆ Dicom Option	Inbuilt
◆ Power Supply	10 KVA	◆ Recommended			Below 25 Cases



Magnetic Resonance Imaging Machine

Matching with the customer's requirements, we are engaged in offering comprehensive range of **Magnetic Resonance Imaging Machine**. We offer MRI Machines to our patrons as per the need of the application. The offered MRI Machines are known for their long functional life and accurate results. Moreover, our MRI Machines are reliable too.

Channel	2 Ch	Tesla	0.3 Tesla	Gap	43 cm V
Magnet	Permanent	Monitor type	CRT		
Magnet Type	Open Type	Dicom Option	Optional		
Gradient Strength	8mt/m	Power Supply	10 KVA		
Slew Rate	11mt/m/Sec	Recommended	Below 15 Cases		



Magnetic Resonance Imaging

Being pioneers in the industry, we are engaged in offering extensive array of **Magnetic Resonance Imaging**. Our products control energy deposition in patients. We offer MRI Machines to our clients as per their needs and in accurate dimensions. Besides these, the offered products are technically advanced and available at industry leading prices.

Channel	4 Ch	Tesla	0.3 Tesla	Gap	42.5 cm V
Magnet	Permanent	Monitor type	TFT (LCD)		
Magnet Type	Open Type	Dicom Option	Inbuilt		
Gradient Strength	21mt/m	Power Supply	10 KVA		
Slew Rate	55mt/m/Sec	Recommended	Below 25 Cases		



MRI Systems

Owing to years of experience in domain, we are engaged in offering wide array of MRI Systems. The offered products are used for magnetic navigation, improving the sonographic identification and many more purposes. We offer MRI Systems to our clients in recommended specifications and used in various places. Moreover, we offer MRI Systems at affordable rates.

Slew Rate	7mt/m/Sec	Tesla	0.3 Tesla	Channel	2 Ch
Magnet	Permanent	Gap	41 cm V	Magnet Type	Close Type
Monitor type	CRT	Gradient Strength	8mt/m		
Dicom Option	Optional	Power Supply	10 KVA		




CT Scanners Machines

We are offering wide range of **CT Scanners Machines for Operation Theatres**. Our products use special x-ray equipment to make cross-sectional pictures of the body. We provide CT Scan Machines to our clients as per the need of the application and at industry leading prices. Moreover, our CT Scan Machines are reliable too.

Features:

- Technologically advanced
- High performance
- Low maintenance



C Arm X Ray Machines

We are dedicatedly engaged in trading and supplying high quality **C Arm X Ray Machines**. Manufactured using high quality component and advanced technology, these X-ray machines are in tandem with international quality standards. These machines are highly demanded in medical field due to their various applications and excellent performance.



Features:

- High performance
- Low maintenance required
- Tough construction

Mammography Units

Customers can avail **Mammography Units** from us as per the industry standards. We offer Mammography Units to our customers as per their needs and used for diagnosis internal diseases of patients. Customers like our Mammography Units for their hassle free performance and best quality. Moreover, we offer Mammography Units at industry leading prices.

Features:

- High demand
- Effective
- Easy to use



Ultrasound Transducers

Leveraging the vast industrial experience, we are engaged in offering bulk of **Ultrasound Transducers**. These Ultrasound Transducers are extensively used in the manufacturing of Ultrasound machines. Also, these Ultrasound Transducers are manufactured by using advance technology and high grade material. Also, these Ultrasound Transducers are known for their superior quality.

- Longer service life
- Efficient working
- High quality



Magnetic Resonance Imaging Machine for Hospital

We are offering wide range of **Magnetic Resonance Imaging Machine for Hospital**. We offer MRI Machines to our patrons as per the need of the application. The offered MRI Machines are known for their long functional life and accurate results.

Features:

- High quality
- Effectiveness
- High performance



X Ray Machines

We are offering wide range of **X Ray Machines for Medical Laboratories**. Our products are used in hospitals for diagnosis various body parts. The offered products are known for their accuracy and reliability. The demand of these X Ray Machines is highly demanded in medical field due to their various applications and excellent performance.

Features:

- ◆ Effectiveness
- ◆ Premium quality
- ◆ High performance



DEFIBRILLATOR FEATURES

- Large 8.4" colour TFT screen with LED backlight
- Single rotary knob operation
- Direct function keys
- Visual instructions and guidance on screen
- 3 ECG strips memory
- 24 hours trend for all parameters
- Adult and pediatric application
- Lightweight & portable



Dialysis Machines

Quality treatment to suit your budget quality

Almost one in two patients with ESRD dies as a result of cardiovascular disease. That is why Cardioprotective Hemodialysis is a core principle of Fresenius Medical Care, as we work and strive to solve the challenges of modern dialysis. Each step we take is focused on minimizing cardiovascular risks and extending patients' lives. The Next Generation 4008S continues the success story of the well-known 4008 series, by combining best quality HD treatment, proven reliability and operational efficiency.

Proven quality

The latest member of a long established and highly successful 4008 series

- Advanced dialysis fluid circuit with highly-precise volumetrically controlled ultrafiltration
- Consistently reliable hydraulics

Advanced treatment

High-quality HD treatment to protect your patient

- OCM® (Online Clearance Monitor)
- DIASAFE® plus (Dialysis Fluid Filter)
- bibag® (Dry Bicarbonate Concentrate)
- Single needle or SN Click-Clack
- ISO-UF program



PATIENT MONITORS

FEATURES

- 12.1" high resolution color TFT LCD display
- Critical parameters ECG, SpO₂, NIBP, Temp, Resp monitored
- 12 lead ECG waveform display & total 11 waveforms display
- ST analysis, arrhythmia analysis
- Digital SpO₂ even for low perfusion cases
- 72 hours tabular and graphical trends & data storage
- Drug dose and Hemodynamic calculation software
- OXY CRG software
- CNS : Wired
- Certification: CE 0123, FDA



MULTIPARA PATIENT MONITOR

MPM 212

- § 12.1" TFT touch screen cluster-free multi-display formats
- § Intelligent alarm management for effective patient monitoring
- § Non-volatile storage data (offline memory)
- § 96-hour trends & graphical recall, 60 alarm events
- § Networking capability (wired/wireless)
- § Upgrade options for Thermal Recorder, Dual IBP, EtcO₂, CO, AGM, ICG, IOC



Digital X Ray



Ventilator

www.kiransurgicalsystems.com



Delhi : A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092

Gujarat : 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda - 390 010,

Andhra Pradesh : No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad - 500 034

Tamil Nadu : No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore -641 021

Kerala : 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621 ,

Helpline : +91 9751725000, +91 9751625000, + 91 9751425000, + 919751215000, + 91 9751245000 + 91 9846596000

Email : kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com

Website : www.kirangroupofcompanies.com, www.kiransurgicalsystems.com, www.kiran.biz, www.kiranincenerationsystems.com
www.kiransolarsystems.com, www.kiranled.com, www.kiranequipments.com, www.kiranpump.com
www.kiranpanels.com, www.kiranappliances.com, www.kirancompressors.com, www.kiranindia.in