


☐ GENERAL INSTRUMENTS



□ CONTAINERS - S.S.



**Container for 3.5 & 4mm Screws
with Two Instruments Trays**
Code No.: 400.010



Container for 3.5 & 4.0mm Screws
Code No.: 400.025



**Container for 3.5 & 4.0mm Screws,
Large with One Instrument Tray**
Code No.: 400.026



Container for 4.5 & 6.5mm Screws
Code No.: 400.020




**Container with Two Trays
for Instruments, 300mm long**
Code No.: 400.005

S.S.CONTAINERS

Code No.	Description
400.005	Container with Two Trays for Instruments, 300mm long
400.010	Container for 3.5 & 4mm Screws with Two Instruments Trays
400.020	Container for 4.5 & 6.5mm Screws
400.025	Container for 3.5 & 4.0mm Screws
400.026	Container for 3.5 & 4.0mm Screws, Large with One Instrument Tray

□ STERILIZATION CONTAINER - PLASTIC WITH SILICONE MAT



STERILIZATION CONTAINER, PLASTIC WITH SILICONE MAT

Item Code	Description	Width inch (mm)	Length inch (mm)	Height inch (mm)
491.101	Single Tier Tray with 1 Silicone Mat	2.5(64)	6.0(152)	0.75(19)
491.105	Single Tier Tray with 1 Silicone Mat	4.0(102)	7.5(190)	0.75(19)
491.109	Single Tier Tray with 1 Silicone Mat	6.0(152)	10.0(254)	0.75(19)
491.113	Deep Single Tier Tray with 1 Silicone Mat	6.0(152)	10.0(254)	2.25(57)
491.117	Deep Double Tiered Tray with 2 Silicone Mats	6.0(152)	10.0(254)	2.25(57)



BONE TAP
WITH T-HANDLE

Code No.	Dia. mm	
263.010	1.5	Cortex
263.015	2.0	Cortex
263.020	2.7	Cortex
263.025	3.5	Cortex
263.030	3.5	Cancellous
263.038	4.0	Cancellous
263.040	4.5	Cortex
263.045	4.5	Cortex
		Extra Long
263.050	6.5	Cancellous
		Long



BONE TAP
QUICK COUPLING END

Code No.	Dia. mm	
264.020	2.7	Cortex
264.025	3.5	Cortex
264.030	3.5	Cancellous
264.038	4.0	Cancellous
264.040	4.5	Cortex
264.045	4.5	Cortex
		Extra Long
264.050	6.5	Cancellous
		Long

GENERAL INSTRUMENTS

...DRILL BITS, BONE TAPS, REAMERS & CLAMPS

YouTube www.TV.ortho.in/10180

YouTube www.TV.ortho.in/25158

DRILL BIT - QUICK COUPLING END

Dia. mm	Length mm	Codes	
		S.S.	Titanium Nitrate
1.5	60	260.150	TN-260.150
1.5	100	260.151	TN-260.151
2.0	100	260.200	TN-260.200
2.5	115	260.250	TN-260.250
2.5	130	260.251	TN-260.251
2.5	150	260.252	TN-260.252
2.7	150	260.270	TN-260.270
2.7	130	260.271	TN-260.271
2.7	200	260.273	TN-260.273
2.8	165	260.280	TN-260.280
2.8	180	260.282	TN-260.282
3.0	130	260.300	TN-260.300
3.0	150	260.301	TN-260.301
3.0	200	260.303	TN-260.303
3.2	130	260.320	TN-260.320
3.2	150	260.321	TN-260.321
3.2	200	260.325	TN-260.325
3.2	250	260.327	TN-260.327
3.2	300	260.330	TN-260.330
3.5	130	260.350	TN-260.350
3.5	150	260.351	TN-260.351
3.5	200	260.375	TN-260.375
3.5	250	260.376	TN-260.377
3.5	300	260.380	TN-260.380
4.0	130	260.400	TN-260.400
4.0	150	260.401	TN-260.401
4.0	200	260.405	TN-260.405
4.0	250	260.407	TN-260.407
4.0	300	260.410	TN-260.410
4.3	221	260.430	TN-260.430
4.3	210	260.431	TN-260.431
4.3	280	260.432	TN-260.432
4.5	130	260.450	TN-260.450
4.5	150	260.451	TN-260.451
4.5	200	260.475	TN-260.475
4.5	250	260.477	TN-260.477
4.5	300	260.480	TN-260.480
4.0 x 4.5	250	260.501	TN-260.501
5.0	250	260.502	TN-260.502

DRILL BIT - PLAIN SHANK / JACOB CHUCK END

Dia. mm	Length mm	Codes	
		S.S.	Titanium Nitrate
1.1	60	261.190	TN-261.190
1.5	60	261.192	TN-261.192
1.5	100	261.191	TN-261.191
2.0	100	261.620	TN-261.620
2.5	115	261.625	TN-261.625
2.5	130	261.626	TN-261.626
2.5	150	261.621	TN-261.621
2.7	150	261.622	TN-261.622
2.7	130	261.627	TN-261.627
2.7	200	261.628	TN-261.628
2.8	165	261.623	TN-261.623
2.8	180	261.624	TN-261.624
3.0	130	261.630	TN-261.630
3.0	150	261.629	TN-261.629
3.0	200	261.631	TN-261.631
3.2	130	261.632	TN-261.632
3.2	150	261.634	TN-261.634
3.2	200	261.633	TN-261.633
3.2	250	261.638	TN-261.638
3.2	300	261.639	TN-261.639
3.5	130	261.635	TN-261.635
3.5	150	261.641	TN-261.641
3.5	200	261.636	TN-261.636
3.5	250	261.637	TN-261.637
3.5	300	261.642	TN-261.642
4.0	130	261.640	TN-261.640
4.0	150	261.643	TN-261.643
4.0	200	261.644	TN-261.644
4.0	250	261.647	TN-261.647
4.0	300	261.648	TN-261.648
4.3	221	261.653	TN-261.653
4.3	280	261.654	TN-261.654
4.5	130	261.645	TN-261.645
4.5	150	261.649	TN-261.649
4.5	200	261.646	TN-261.646
4.5	250	261.651	TN-261.651
4.5	300	261.652	TN-261.652
5.0	250	261.650	TN-261.650
6.0	250	261.660	TN-261.660

YouTube www.TV.ortho.in/25175



CANNULATED REAMER (T-TYPE)

Code No.	Size
507.96-8	8 mm
507.96-9	9 mm
507.96-10	10 mm
507.96-11	11 mm
507.96-12	12 mm
507.96-13	13 mm
507.96-14	14 mm



LOWMAN BONE CLAMP - STAINLESS STEEL

Code No.	Description
352.010	for Fingers (Mini)
352.020	for Forearms (Extra Small)
352.030	for Humerus (Small)
352.040	for Tibia (Medium)
352.050	for Femur (Large)

YouTube www.TV.ortho.in/4021



PATELLA BONE HOLDING FORCEPS, THREE PRONG
Code No.: 343.010



PATELLA BONE HOLDING FORCEPS, FOUR PRONG
Code No.: 343.020



BONE AND PLATE HOLDING FORCEPS
Code No.: 348.010



MODIFIED PATELLA BONE HOLDING FORCEPS FOUR PRONG
Code No.: 343.024



MALLEOLAR FORCEPS WITH 2 POINTS WHICH MEET AT 15° ANGLE 210MM
Code No.: 347.000

YouTube www.TV.ortho.in/4480



MODIFIED PATELLA BONE HOLDING FORCEPS THREE PRONG
Code No.: 343.022



BURN'S BONE HOLDING FORCEPS HEAVY DUTY
Code No.: 350.000



FIBULAR FORCEPS - STRAIGHT
Code No.: 345.010



LANE BONE HOLDING FORCEPS WITH RATCHET

Code No.	Length mm
349.010	200
349.020	250
349.030	300



LANE BONE HOLDING FORCEPS WITHOUT RATCHET LOCK

Code No.	Length mm
349.050	200
349.060	250
349.070	300



**REDUCTION FORCEPS
SERRATED, SPEED LOCK**

Code No.	Length mm
341.025	140
341.030	170
341.040	200
341.050	240
341.055	280

[YouTube www.TV.ortho.in/25183](https://www.youtube.com/watch?v=www.TV.ortho.in/25183)



**REDUCTION FORCEPS,
POINTED, SPEED LOCK**

Code No.	Length mm
342.035	140
342.040	170
342.050	200
342.060	240
342.065	280

[YouTube www.TV.ortho.in/25210](https://www.youtube.com/watch?v=www.TV.ortho.in/25210)



**REDUCTION FORCEPS
SERRATED, RATCHET LOCK**

Code No.	Length mm
341.010	140
341.020	170
341.022	200



**REDUCTION FORCEPS
POINTED, RATCHET LOCK**

Code No.	Length mm
342.010	140
342.020	170
342.030	200

**SELF CENTERING BONE
HOLDING FORCEPS**

Code No.	Length mm
340.010	150
340.020	190
340.030	240
340.040	260
340.050	280



PELVIC REDUCTION FORCEPS, S.S.

- 344.010 Pelvic Reduction Forceps, Oblique, Small, Length-190mm, with Pointed Ball Tips, Speed Lock
- 344.020 Pelvic Reduction Forceps, Oblique, Small, Length-240mm, with Pointed Ball Tips, Speed Lock
- 344.030 Pelvic Reduction Forceps, Length-400mm with Pointed ball Tips, Speed Lock
- 344.040 Pelvic Reduction Forceps, Small for Screws Length-200mm, use with 3.5 and 4.5mm Cortex Screws
- 344.050 Pelvic Reduction Forceps, Medium for Screws Length-250mm, use with 3.5 and 4.5mm Cortex Screws
- 344.060 Pelvic Reduction Forceps, Large for Screws Length-330mm, Speed Lock, Fixation with Two 4.5mm Cortex Screws
- 344.070 Pelvic Reduction Forceps, with Angled Pointed Ball Tips Length-250mm, Speed Lock
- 344.080 Pelvic Reduction Forceps, Long, with Pointed Ball Tips Length-400mm, Speed Lock
- 344.090 Pelvic Reduction Forceps, Long, with Three Pointed Ball Tips Length-400mm, Speed Lock

FIBULAR FORCEPS, S.S.

- 345.010 210mm - Straight
- 345.020 210mm - Curved

HOLDING FORCEPS FOR TIBIAL EDGE FRAGMENT, S.S.

- 346.000 210mm

- 351.000 Doyen's Rib Shearer, S.S.



**BONE NIBBLER
(DOUBLE ACTION) - CURVED**

Code No.	Length mm
356.010	190
356.020	225
356.030	270

YouTube www.TV.ortho.in/25233



**BONE NIBBLER
(DOUBLE ACTION) - STRAIGHT**

Code No.	Length mm
355.010	190
355.020	225
355.030	270



**BONE NIBBLER
(DOUBLE ACTION) - ANGULAR**

Code No.	Length mm
357.010	190
357.020	225
357.030	270

YouTube www.TV.ortho.in/25238



**TUDOR EDWARD
BONE CUTTING FORCEPS**

Code No.: 358.000

YouTube www.TV.ortho.in/5588



**BONE CUTTING FORCEPS
(DOUBLE ACTION) - STRAIGHT**

Code No.	Length mm
353.010	190
353.020	225
353.030	270



**BONE CUTTING FORCEPS
(DOUBLE ACTION) - CURVED**

Code No.	Length mm
354.010	190
354.020	225
354.030	270

AS023 BONE HOLDING FORCEPS SET, S.S.

Set Consisting of:

		Units
340.030	Self Centering Bone Holding Forceps - 240mm	1
340.040	Self Centering Bone Holding Forceps - 260mm	2
340.050	Self Centering Bone Holding Forceps - 280mm	2
341.030	Reduction Forceps Serrate Jaw Speed Lock - 170mm.	1
341.050	Reduction Forceps Serrate Jaw Speed Lock - 240mm.	1
342.040	Reduction Forceps Pointed Speed Lock - 170mm	1
342.050	Reduction Forceps Pointed Speed Lock - 200mm	2
903.023	Empty Container	

GENERAL INSTRUMENTS

...BONE FILE, GAUGE, OSTEOTOMES & BONE CURETTE



**BONE FILE, S.S.
WITH FIBRE HANDLE FLAT**

Code No.	Size mm
361.020	100
361.030	150

[You Tube www.TV.ortho.in/25249](https://www.youtube.com/watch?v=www.TV.ortho.in/25249)



BONE FILE, HALF ROUND, POINTED
Code No.: 361.010



**OSTEOTOMES, S.S.
WITH FIBRE HANDLE STRAIGHT**

Code No.	Size mm
363.005	5
363.010	10
363.015	15
363.020	20
363.025	25
363.030	30

[You Tube www.TV.ortho.in/25273](https://www.youtube.com/watch?v=www.TV.ortho.in/25273)



**GOUGE, S.S.
WITH FIBRE HANDLE STRAIGHT**

Code No.	Size mm
365.005	5
365.010	10
365.015	15
365.020	20
365.025	25
365.030	30

[You Tube www.TV.ortho.in/25275](https://www.youtube.com/watch?v=www.TV.ortho.in/25275)



**BONE CURETTE
DOUBLE ENDED**
Code No.: 367.010

[You Tube www.TV.ortho.in/4602](https://www.youtube.com/watch?v=www.TV.ortho.in/4602)



**CHISELS, S.S.
WITH FIBRE HANDLE**

Code No.	Size mm
362.005	5
362.010	10
362.015	15
362.020	20
362.025	25
362.030	30

[You Tube www.TV.ortho.in/25264](https://www.youtube.com/watch?v=www.TV.ortho.in/25264)



**OSTEOTOMES, S.S.
WITH FIBRE HANDLE CURVED**

Code No.	Size mm
364.005	5
364.010	10
364.015	15
364.020	20
364.025	25
364.030	30

[You Tube www.TV.ortho.in/25274](https://www.youtube.com/watch?v=www.TV.ortho.in/25274)



**GOUGE, S.S.
WITH FIBRE HANDLE CURVED**

Code No.	Size mm
366.005	5
366.010	10
366.015	15
366.020	20
366.025	25
366.030	30

[You Tube www.TV.ortho.in/25276](https://www.youtube.com/watch?v=www.TV.ortho.in/25276)



**BONE CURETTE
WITH FIBRE HANDLE**
Code No.: 367.020



BONE HAMMER - S.S.

Code No.	Weight gms
368.010	200
368.020	350
368.030	500
368.040	700

YouTube www.TV.ortho.in/14560



BONE HOOK

Code No.	Size
369.010	Small
369.020	Medium
369.030	Large

YouTube www.TV.ortho.in/25283



FINGER SAW

Code No.: 360.000



SHARP HOOK

Code No.: 371.000



AMPUTATION SAW

Code No.: 359.000



BONE HAMMER - NYLON FACED

Code No.: 368.050

YouTube www.TV.ortho.in/7959



SKIN HOOK - 1 PRONG

Code No.: 370.010



SKIN HOOK - 2 PRONG

Code No.: 370.020



SKIN HOOK - 3 PRONG

Code No.: 370.030



LANE'S BONE LEVER SERRATED

Code No.: 372.010



**LANE'S BONE LEVER, PLAIN, SMALL
STAINLESS STEEL**

Code No.: 372.020



**LANE'S BONE LEVER, PLAIN, LARGE
STAINLESS STEEL**

Code No.: 372.030

AS024 GENERAL INSTRUMENTS SET, S.S.

Set Consisting of:

	Units
369.010 Bone Hook - Small	1
369.020 Bone Hook - Medium	1
379.025 Hohmann Retractor 8mm Wide, Short Narrow Tip, Length 240mm	2
379.030 Hohmann Retractor 18mm Wide, Short Narrow Tip, Length 240mm	2
379.045 Hohmann Retractor 24mm Wide, with Long wide Tip, Length 240mm	1
381.012 Periosteal Elevator Straight with Fibre Handle - 12mm	2
368.030 Hammer - 500gm	1
368.050 Hammer Nylon Faced	1
366.015 Bone Gauge Curved with Fibre Handle, 15mm	1
903.025 Empty Container	



**BONE SPREADER
STAINLESS STEEL**

Code No.	Length mm
383.014	140
383.021	210
383.027	270



www.TV.ortho.in/25294

**LAGENBACK RETRACTOR
STAINLESS STEEL**

Code No.	Size
374.010	Small
374.015	Medium
374.020	Large

www.TV.ortho.in/14130



**BRISTOW'S RETRACTOR
STAINLESS STEEL**

Code No.: 373.000



**HOHMANN RETRACTOR
240MM**



**HOHMANN RETRACTOR WIDE TIP
240MM**



**HOHMANN RETRACTOR
160MM**

Code No.	Description
379.010	Hohmann Retractor 6mm Wide, Short Narrow Point, Length 160mm
379.015	Hohmann Retractor 8mm Wide, Short Narrow Point, Length 160mm
379.020	Hohmann Retractor 15mm Wide, Length 160mm
379.025	Hohmann Retractor 8mm Wide, Short Narrow Tip, Length 240mm
379.030	Hohmann Retractor 18mm Wide, Short Narrow Tip, Length 240mm
379.035	Hohmann Retractor 18mm Wide, Long Narrow Tip, Length 240mm
379.040	Hohmann Retractor 22mm Wide, with Wide Tip, Length 240mm
379.045	Hohmann Retractor 24mm Wide, with Long wide Tip, Length 240mm
379.050	Hohmann Retractor 43mm Wide, Narrow Tip, Length 240mm
379.055	Hohmann Retractor 70mm Wide, Narrow Tip, Length 240mm



VOLKMAN'S RETRACTOR

Code No.	No. of Prong
375.010	2
375.020	3
375.030	4
375.040	5

www.TV.ortho.in/9835



RIGHT ANGLE RETRACTOR

Code No.: 376.000



DOUBLE HOOK RETRACTOR

Code No.: 377.000

www.TV.ortho.in/15121

SELF RETAINING RETRACTOR

Code No.	Size
380.010	Small
380.020	Medium
380.030	Large
380.040	Extra Large

YouTube www.TV.ortho.in/27076



SKIN RETRACTOR
Code No.: 377.050



**PERIOSTEAL ELEVATOR - STRAIGHT
STAINLESS STEEL**

Code No.	Length mm
381.106	6
381.108	8
381.110	10
381.112	12
381.114	14
381.116	16



**PERIOSTEAL ELEVATOR - CURVED
STAINLESS STEEL**

Code No.	Length mm
382.106	6
382.108	8
382.110	10
382.112	12
382.114	14
382.116	16



**PERIOSTEAL ELEVATORS- STRAIGHT
WITH FIBRE HANDLE**

Code No.	Length mm
381.006	6
381.008	8
381.010	10
381.012	12
381.014	14
381.016	16

YouTube www.TV.ortho.in/25315



**PERIOSTEAL ELEVATORS - CURVED
WITH FIBRE HANDLE**

Code No.	Length mm
382.006	6
382.008	8
382.010	10
382.012	12
382.014	14
382.016	16

YouTube www.TV.ortho.in/25313

**ELECTRIC PLASTER CUTTING SAW
(220V 50Hz AC)**

Code No.	Description
399.005	Electric Plaster Cutting Saw with Two Blades
399.006	Spare Blade - Small
399.007	Spare Blade - Large





PLASTER SAW - BERGMAN
Code No.: 388.015



PLASTER SAW - HEAVY HANDLE
Code No.: 388.005



PLASTER SAW - ENGEL
Code No.: 388.010



PLASTER SPREADER
Code No.: 388.020



PLASTER SHEAR - GAY'S
Code No.: 388.035



PLASTER SHEAR - BOHLER'S
Code No.: 388.025



PLASTER SHEAR - STILLE'S
Code No.: 388.030



PLASTER OPENER - DAW'S
Code No.: 388.050



PLASTER BENDER
Code No.: 388.045



BANDAGE CUTTING SCISSOR
Code No.: 388.040

□ KNEE & SKULL SURGERY INSTRUMENTS

INSTRUMENTS FOR KNEE SURGERY, S.S.

Code No.	Description
393.005	Blount Knee Retractor
393.010	Martin's Cartilage Clamp
393.015	Maniscus Forceps
393.020	Smilie's Knee Joint Retractor (Set of Six)
393.025	Smillie's Maniscus Knives (Set of Three)
393.030	Charnley's Compression Clamp
393.040	Mullar Compression Clamp



**SMILIE'S KNEE JOINT RETRACTOR
(SET OF SIX)**
Code No.: 393.020




MULLAR COMPRESSION CLAMP
Code No.: 393.040



CHARNLEY'S COMPRESSION CLAMP
Code No.: 393.030

INSTRUMENTS FOR SKULL SURGERY, S.S.

Code No.	Description
394.005	Gigli Saw Wire - Ordinary
394.010	Gigli Saw Wire - Fine (Imported)
394.015	Gigli Saw Handle (In Pair)
394.020	Gigli Saw Introducer
394.025	Crutch Field Tong with 2 Burrs-Small
394.030	Crutch Field Tong with Two Burrs - Medium
394.035	Crutch Field Tong with Two Burrs - Large
394.040	Spare Burrs for Crutch Field Tong
394.045	Gardener - Wells Traction Tongs



**SMILIE'S MANISCUS KNIVES
(SET OF THREE)**
Code No.: 393.025



**GIGLI SAW HANDLE
(IN PAIR)**
Code No.: 394.015
www.TV.ortho.in/14480



**CRUTCH FIELD TONG
WITH TWO BURRS - MEDIUM**
Code No.: 394.030

YouTube www.TV.ortho.in/14478



GIGLI SAW WIRE - ORDINARY
Code No.: 394.005



GIGLI SAW WIRE - FINE
Code No.: 394.010



DISTRACTOR FOR TIBIA & FEMUR - LARGE
Code No.: 396.015



**DISTRACTOR FOR TIBIA & FEMUR
EXTRA LARGE**
Code No.: 396.016



**PAEDIATRIC TIBIAL & FEMUR
REDUCTION DEVICE (F-TOOL SMALL)**
Code No.: 398.042



SKIN GRAFTING KNIFE, EXCEL
Code No.: 395.020



STERILIZATION TUBES

Code No.	Length mm
401.002	150
401.004	225
401.006	300
401.008	400
401.010	500

Profile:
Stainless Steel Tubes for Sterilization of
K. Wires, Nails, Pins & Drill Bits



**TIBIAL & FEMUR REDUCTION DEVICE
(F-TOOL LARGE)**
Code No.: 398.053

ESMARCH RUBBER BANDAGE

Code No.	Description
397.005	Esmarch Rubber Bandage - 75mm
397.010	Esmarch Rubber Bandage - 100mm
397.015	Esmarch Rubber Bandage- 150mm

□ PNEUMATIC TOURNIQUET SET - MANUAL


YouTube www.TV.ortho.in/14937

PNEUMATIC TOURNIQUET - MANUAL

Code No.	Description
398.000	Pneumatic Tourniquet, Manual, alongwith Handpump & 3 Cuffs
398.063	Pneumatic Tourniquet Set, Electronic (220V 50Hz AC) Digital Display with 3 Hrs. Batterybackup.
398.066	Pneumatic Tourniquet Set (Dual), Electronic (220V 50Hz AC) Digital Display with 3 Hrs. Batterybackup.
398.070	Rubber Cuffs, Washable and Easy Fitting, Set of 5 Pcs. - Pediatric, Small, Medium, Large & Extra Large.
398.078	Silicone Cuffs, Washable, Easy Fitting & Autoclavable Set of 5 Pcs. - Pediatric, Small, Medium, Large & Extra Large.



**PNEUMATIC TOURNIQUET
ALONGWITH HANDPUMP & 3 CUFFS**
Code No.: 398.000



SILICONE CUFFS, SET OF 5 PCS.
Code No.: 398.078



RUBBER CUFFS, SET OF 5 PCS.
Code No.: 398.070

□ PNEUMATIC TOURNIQUET SET - ELECTRIC



PNEUMATIC TOURNIQUET SET - ELECTRONIC (220V 50HZ AC),
DIGITAL DISPLAY WITH 3 HRS. BATTERY BACKUP
Code No.: 398.063



PNEUMATIC TOURNIQUET SET (DUAL) - ELECTRONIC (220V 50HZ AC),
DIGITAL DISPLAY WITH 3 HRS. BATTERYBACKUP
Code No.: 398.066