



## □ HAND & FOOT IMPLANTS






**CORTEX SCREW 1.5MM**

| Length in mm | Code S.S. | Code Titanium |
|--------------|-----------|---------------|
| 6            | 118.006   | Ti-118.006    |
| 7            | 118.007   | Ti-118.007    |
| 8            | 118.008   | Ti-118.008    |
| 9            | 118.009   | Ti-118.009    |
| 10           | 118.010   | Ti-118.010    |
| 11           | 118.011   | Ti-118.011    |
| 12           | 118.012   | Ti-118.012    |
| 14           | 118.014   | Ti-118.014    |
| 16           | 118.016   | Ti-118.016    |
| 18           | 118.018   | Ti-118.018    |
| 20           | 118.020   | Ti-118.020    |

Drill Bit: 1.1mm

[YouTube www.TV.ortho.in/24564](https://www.youtube.com/watch?v=www.TV.ortho.in/24564)



**CORTEX SCREW 2.0MM**

| Length in mm | Code S.S. | Code Titanium |
|--------------|-----------|---------------|
| 6            | 118.106   | Ti-118.106    |
| 7            | 118.107   | Ti-118.107    |
| 8            | 118.108   | Ti-118.108    |
| 9            | 118.109   | Ti-118.109    |
| 10           | 118.110   | Ti-118.110    |
| 11           | 118.111   | Ti-118.111    |
| 12           | 118.112   | Ti-118.112    |
| 14           | 118.114   | Ti-118.114    |
| 16           | 118.116   | Ti-118.116    |
| 18           | 118.118   | Ti-118.118    |
| 20           | 118.120   | Ti-118.120    |

Drill Bit: 1.5mm

[YouTube www.TV.ortho.in/156](https://www.youtube.com/watch?v=www.TV.ortho.in/156)



**CORTEX SCREW 1.5MM SELF TAPPING**

| Length in mm | Code S.S. | Code Titanium |
|--------------|-----------|---------------|
| 6            | 118.206   | Ti-118.206    |
| 7            | 118.207   | Ti-118.207    |
| 8            | 118.208   | Ti-118.208    |
| 9            | 118.209   | Ti-118.209    |
| 10           | 118.210   | Ti-118.210    |
| 11           | 118.211   | Ti-118.211    |
| 12           | 118.212   | Ti-118.212    |
| 14           | 118.214   | Ti-118.214    |
| 16           | 118.216   | Ti-118.216    |
| 18           | 118.218   | Ti-118.218    |
| 20           | 118.220   | Ti-118.220    |

Drill Bit: 1.1mm




**CORTEX SCREW 2.0MM SELF TAPPING**

| Length in mm | Code S.S. | Code Titanium |
|--------------|-----------|---------------|
| 6            | 118.306   | Ti-118.306    |
| 7            | 118.307   | Ti-118.307    |
| 8            | 118.308   | Ti-118.308    |
| 9            | 118.309   | Ti-118.309    |
| 10           | 118.310   | Ti-118.310    |
| 11           | 118.311   | Ti-118.311    |
| 12           | 118.312   | Ti-118.312    |
| 14           | 118.314   | Ti-118.314    |
| 16           | 118.316   | Ti-118.316    |
| 18           | 118.318   | Ti-118.318    |
| 20           | 118.320   | Ti-118.320    |

Drill Bit: 1.5mm


YouTube [www.TV.ortho.in/24582](http://www.TV.ortho.in/24582)



**CORTEZ SCREW 2.4MM**

| Length in mm | Code S.S. | Code Titanium |
|--------------|-----------|---------------|
| 6            | 100.006   | Ti-100.006    |
| 8            | 100.008   | Ti-100.008    |
| 10           | 100.010   | Ti-100.010    |
| 12           | 100.012   | Ti-100.012    |
| 14           | 100.014   | Ti-100.014    |
| 16           | 100.016   | Ti-100.016    |
| 18           | 100.018   | Ti-100.018    |
| 20           | 100.020   | Ti-100.020    |
| 22           | 100.022   | Ti-100.022    |
| 24           | 100.024   | Ti-100.024    |
| 26           | 100.026   | Ti-100.026    |
| 28           | 100.028   | Ti-100.028    |
| 30           | 100.030   | Ti-100.030    |

Drill Bit: 1.8mm




**CORTEZ SCREW 2.7MM**

| Length in mm | Code S.S. | Code Titanium |
|--------------|-----------|---------------|
| 6            | 101.006   | Ti-101.006    |
| 8            | 101.008   | Ti-101.008    |
| 10           | 101.010   | Ti-101.010    |
| 12           | 101.012   | Ti-101.012    |
| 14           | 101.014   | Ti-101.014    |
| 16           | 101.016   | Ti-101.016    |
| 18           | 101.018   | Ti-101.018    |
| 20           | 101.020   | Ti-101.020    |
| 22           | 101.022   | Ti-101.022    |
| 24           | 101.024   | Ti-101.024    |
| 26           | 101.026   | Ti-101.026    |
| 28           | 101.028   | Ti-101.028    |
| 30           | 101.030   | Ti-101.030    |

Drill Bit: 2mm


Length 32 to 55mm available on request.



**CORTEZ SCREW 2.4MM SELF TAPPING**

| Length in mm | Code S.S. | Code Titanium |
|--------------|-----------|---------------|
| 6            | 100.206   | Ti-100.206    |
| 8            | 100.208   | Ti-100.208    |
| 10           | 100.210   | Ti-100.210    |
| 12           | 100.212   | Ti-100.212    |
| 14           | 100.214   | Ti-100.214    |
| 16           | 100.216   | Ti-100.216    |
| 18           | 100.218   | Ti-100.218    |
| 20           | 100.220   | Ti-100.220    |
| 22           | 100.222   | Ti-100.222    |
| 24           | 100.224   | Ti-100.224    |
| 26           | 100.226   | Ti-100.226    |
| 28           | 100.228   | Ti-100.228    |
| 30           | 100.230   | Ti-100.230    |

Drill Bit: 1.8mm



**CORTEZ SCREW 2.7MM SELF TAPPING**

| Length in mm | Code S.S. | Code Titanium |
|--------------|-----------|---------------|
| 6            | 101.206   | Ti-101.206    |
| 8            | 101.208   | Ti-101.208    |
| 10           | 101.210   | Ti-101.210    |
| 12           | 101.212   | Ti-101.212    |
| 14           | 101.214   | Ti-101.214    |
| 16           | 101.216   | Ti-101.216    |
| 18           | 101.218   | Ti-101.218    |
| 20           | 101.220   | Ti-101.220    |
| 22           | 101.222   | Ti-101.222    |
| 24           | 101.224   | Ti-101.224    |
| 26           | 101.226   | Ti-101.226    |
| 28           | 101.228   | Ti-101.228    |
| 30           | 101.230   | Ti-101.230    |

Drill Bit: 2mm

Length 32 to 55mm available on request.

YouTube [www.TV.ortho.in/24584](http://www.TV.ortho.in/24584)

YouTube [www.TV.ortho.in/157](http://www.TV.ortho.in/157)

YouTube [www.TV.ortho.in/24602](http://www.TV.ortho.in/24602)




**STRAIGHT PLATE 2.0**

| No. of Holes | Length in mm | S.S.    | Code | Titanium   |
|--------------|--------------|---------|------|------------|
| 2            | 11           | 121.102 |      | Ti-121.102 |
| 3            | 17           | 121.103 |      | Ti-121.103 |
| 4            | 23           | 121.104 |      | Ti-121.104 |
| 5            | 29           | 121.105 |      | Ti-121.105 |
| 6            | 35           | 121.106 |      | Ti-121.106 |
| 7            | 41           | 121.107 |      | Ti-121.107 |
| 8            | 47           | 121.108 |      | Ti-121.108 |



**DYNAMIC COMPRESSION PLATE 2.0**


| No. of Holes | Length in mm | S.S.    | Code | Titanium   |
|--------------|--------------|---------|------|------------|
| 4            | 22           | 121.554 |      | Ti-121.554 |
| 5            | 27           | 121.555 |      | Ti-121.555 |
| 6            | 32           | 121.556 |      | Ti-121.556 |
| 7            | 37           | 121.557 |      | Ti-121.557 |
| 8            | 42           | 121.558 |      | Ti-121.558 |



**T PLATE 2.0**

| Length in mm | Code No. S.S. | Length in mm | Code No. TITANIUM |
|--------------|---------------|--------------|-------------------|
| 18           | 123.100       | 18           | Ti-123.100        |

Profile: Shaft with two holes & head with two holes.




**T PLATE 2.7**

| Length in mm | Code No. S.S. | Length in mm | Code No. TITANIUM |
|--------------|---------------|--------------|-------------------|
| 32           | 134.100       | 32           | Ti-134.100        |

Profile: Shaft with three holes & head with two holes.


YouTube [www.TV.ortho.in/112](http://www.TV.ortho.in/112)

YouTube [www.TV.ortho.in/5727](http://www.TV.ortho.in/5727)



**L PLATE 2.7**

| Type          | Length in mm | S.S.    | Code | TITANIUM   |
|---------------|--------------|---------|------|------------|
| Left          | 32           | 133.101 |      | Ti-133.101 |
| Right         | 32           | 133.102 |      | Ti-133.102 |
| Oblique Left  | 34           | 133.103 |      | Ti-133.103 |
| Oblique Right | 34           | 133.104 |      | Ti-133.104 |



**L PLATE 2.0**

| Type          | Length in mm | S.S.    | Code | TITANIUM   |
|---------------|--------------|---------|------|------------|
| Left          | 17           | 122.101 |      | Ti-122.101 |
| Right         | 17           | 122.102 |      | Ti-122.102 |
| Oblique Left  | 19           | 122.103 |      | Ti-122.103 |
| Oblique Right | 19           | 122.104 |      | Ti-122.104 |

Profile: Shaft with two holes & head with two holes.



YouTube [www.TV.ortho.in/68](http://www.TV.ortho.in/68)


#### CONDYLAR PLATE



| Screw in mm | Type      | Length in mm | Code S.S. |
|-------------|-----------|--------------|-----------|
| 1.5         | Pin Left  | 36           | 74.01     |
| 1.5         | Pin Right | 36           | 74.02     |
| 2.0         | Pin Left  | 39           | 74.03     |
| 2.0         | Pin Right | 39           | 74.04     |
| 2.7         | Pin Left  | 57           | 74.05     |
| 2.7         | Pin Right | 57           | 74.06     |

| Screw in mm | Type      | Length in mm | Code S.S. |
|-------------|-----------|--------------|-----------|
| 1.5         | Pin Left  | 36           | Ti-74.01  |
| 1.5         | Pin Right | 36           | Ti-74.02  |
| 2.0         | Pin Left  | 39           | Ti-74.03  |
| 2.0         | Pin Right | 39           | Ti-74.04  |
| 2.7         | Pin Left  | 57           | Ti-74.05  |
| 2.7         | Pin Right | 57           | Ti-74.06  |

YouTube [www.TV.ortho.in/24606](http://www.TV.ortho.in/24606)



#### DYNAMIC COMPRESSION PLATE 2.7

| No. of Holes | Length in mm | Code S.S. |
|--------------|--------------|-----------|
| 2            | 20           | 132.102   |
| 3            | 28           | 132.103   |
| 4            | 36           | 132.104   |
| 5            | 44           | 132.105   |
| 6            | 52           | 132.106   |
| 7            | 60           | 132.107   |
| 8            | 68           | 132.108   |
| 9            | 76           | 132.109   |
| 10           | 84           | 132.110   |
| 12           | 100          | 132.112   |

Thickness (upto 6 Holes): 2mm  
(7 to 12 Holes): 2.5mm

| No. of Holes | Length in mm | Code Titanium |
|--------------|--------------|---------------|
| 2            | 20           | Ti-132.102    |
| 3            | 28           | Ti-132.103    |
| 4            | 36           | Ti-132.104    |
| 5            | 44           | Ti-132.105    |
| 6            | 52           | Ti-132.106    |
| 7            | 60           | Ti-132.107    |
| 8            | 68           | Ti-132.108    |
| 9            | 76           | Ti-132.109    |
| 10           | 84           | Ti-132.110    |
| 12           | 100          | Ti-132.112    |

Thickness (upto 6 Holes): 2mm  
(7 to 12 Holes): 2.5mm

YouTube [www.TV.ortho.in/24608](http://www.TV.ortho.in/24608)



#### QUARTER TUBULAR PLATE 2.7

| No. of Holes | Length in mm | Code S.S. |
|--------------|--------------|-----------|
| 3            | 23           | 131.103   |
| 4            | 31           | 131.104   |
| 5            | 39           | 131.105   |
| 6            | 47           | 131.106   |
| 7            | 55           | 131.107   |
| 8            | 63           | 131.108   |

| No. of Holes | Length in mm | Code Titanium |
|--------------|--------------|---------------|
| 3            | 23           | Ti-131.103    |
| 4            | 31           | Ti-131.104    |
| 5            | 39           | Ti-131.105    |
| 6            | 47           | Ti-131.106    |
| 7            | 55           | Ti-131.107    |
| 8            | 63           | Ti-131.108    |

YouTube [www.TV.ortho.in/158](http://www.TV.ortho.in/158)




RECONSTRUCTION PLATE 2.7 STRAIGHT

| No. of Holes | Length in mm | Code S.S. |
|--------------|--------------|-----------|
| 6            | 48           | 135.106   |
| 8            | 64           | 135.108   |
| 10           | 80           | 135.110   |
| 12           | 96           | 135.112   |
| 14           | 112          | 135.114   |
| 16           | 128          | 135.116   |
| 18           | 144          | 135.118   |
| 20           | 160          | 135.120   |
| 22           | 176          | 135.122   |
| 24           | 192          | 135.124   |

| No. of Holes | Length in mm | Code TITANIUM |
|--------------|--------------|---------------|
| 6            | 48           | Ti-135.106    |
| 8            | 64           | Ti-135.108    |
| 10           | 80           | Ti-135.110    |
| 12           | 96           | Ti-135.112    |
| 14           | 112          | Ti-135.114    |
| 16           | 128          | Ti-135.116    |
| 18           | 144          | Ti-135.118    |
| 20           | 160          | Ti-135.120    |
| 22           | 176          | Ti-135.122    |
| 24           | 192          | Ti-135.124    |

## □ MINI FRAGMENTS IMPLANTS & INSTRUMENTS SET, S.S.



Code No.: AS510

YouTube [www.TV.ortho.in/24618](http://www.TV.ortho.in/24618)

## AS510 MINI FRAGMENTS IMPLANTS & INSTRUMENTS SET, S.S.

| Set Consisting of |                                        |                        |   | Units | Set Consisting of |                                                   |                                   |  | Units |
|-------------------|----------------------------------------|------------------------|---|-------|-------------------|---------------------------------------------------|-----------------------------------|--|-------|
| 118.006           | Cortex Screw,                          | 1.5mm, SS, Length 6mm  | 4 |       | 185.101           | Kirschner Wire with Trocar Tip,                   | 0.8mm Dia, Length 70mm, SS        |  | 10    |
| 118.007           | Cortex Screw,                          | 1.5mm, SS, Length 7mm  | 4 |       | 185.116           | Kirschner Wire with Trocar Tip,                   | 1.0mm Dia, Length 150mm, SS       |  | 10    |
| 118.008           | Cortex Screw,                          | 1.5mm, SS, Length 8mm  | 4 |       | 185.127           | Kirschner Wire with Trocar Tip,                   | 1.2mm Dia, Length 150mm, SS       |  | 10    |
| 118.009           | Cortex Screw,                          | 1.5mm, SS, Length 9mm  | 4 |       | 185.146           | Kirschner Wire with Trocar Tip,                   | 1.6mm Dia, Length 150mm, SS       |  | 10    |
| 118.010           | Cortex Screw,                          | 1.5mm, SS, Length 10mm | 4 |       | 311.002           | Mini Hexagonal Screws Driver                      | with Sleeve (1.5mm), SS           |  | 1     |
| 118.011           | Cortex Screw,                          | 1.5mm, SS, Length 11mm | 4 |       | 311.003           | Mini Drill Sleeve, SS                             |                                   |  | 1     |
| 118.012           | Cortex Screw,                          | 1.5mm, SS, Length 12mm | 4 |       | 311.004           | Mini Tap Sleeve, SS                               |                                   |  | 1     |
| 118.014           | Cortex Screw,                          | 1.5mm, SS, Length 14mm | 4 |       | 311.005           | Mini Depth Gouge, SS                              |                                   |  | 1     |
| 118.016           | Cortex Screw,                          | 1.5mm, SS, Length 16mm | 4 |       | 311.006           | Mini Self Retaining Retractor for Fingers, SS     |                                   |  | 1     |
| 118.106           | Cortex Screw,                          | 2.0mm, SS, Length 6mm  | 4 |       | 311.007           | Mini Self Retaining Retractor for Metacarpals, SS |                                   |  | 1     |
| 118.108           | Cortex Screw,                          | 2.0mm, SS, Length 8mm  | 4 |       | 311.008           | Mini Self Retaining Retractor                     | for Metacarpals ( with Hinge), SS |  | 1     |
| 118.110           | Cortex Screw,                          | 2.0mm, SS, Length 10mm | 4 |       | 311.009           | Mini Neutral & Loaded Drill Guide, SS             |                                   |  | 1     |
| 118.112           | Cortex Screw,                          | 2.0mm, SS, Length 12mm | 4 |       | 311.010           | Mini Plate Bender ( Pair), SS                     |                                   |  | 1     |
| 118.114           | Cortex Screw,                          | 2.0mm, SS, Length 14mm | 4 |       | 311.011           | Mini Screws holding Forceps                       | (for 1.5 & 2.0mm Screws), SS      |  | 1     |
| 118.116           | Cortex Screw,                          | 2.0mm, SS, Length 16mm | 4 |       | 311.012           | Hollow Mill for 1.5mm & 2.0mm Screws, SS          |                                   |  | 1     |
| 118.118           | Cortex Screw,                          | 2.0mm, SS, Length 18mm | 4 |       | 311.013           | Mini Plate Holding Forceps, SS                    |                                   |  | 1     |
| 118.120           | Cortex Screw,                          | 2.0mm, SS, Length 20mm | 4 |       | 311.014           | Mini Plate & Bone Holding Forceps, SS             |                                   |  | 1     |
| 101.006           | Cortex Screw,                          | 2.7mm, SS, Length 6mm  | 5 |       | 311.015           | Mini Reduction Forceps, SS                        |                                   |  | 1     |
| 101.008           | Cortex Screw,                          | 2.7mm, SS, Length 8mm  | 5 |       | 311.016           | Mini Reduction Forceps ( Pointed), SS             |                                   |  | 1     |
| 101.010           | Cortex Screw,                          | 2.7mm, SS, Length 10mm | 5 |       | 311.017           | Mini Periosteal Elevator, SS                      |                                   |  | 1     |
| 101.012           | Cortex Screw,                          | 2.7mm, SS, Length 12mm | 5 |       | 311.018           | Mini Chisels, Size: 4mm, SS                       |                                   |  | 1     |
| 101.014           | Cortex Screw,                          | 2.7mm, SS, Length 14mm | 5 |       | 311.019           | Mini Chisels, Size: 6mm, SS                       |                                   |  | 1     |
| 101.016           | Cortex Screw,                          | 2.7mm, SS, Length 16mm | 5 |       | 311.020           | Mini Chisels, Size: 8mm, SS                       |                                   |  | 1     |
| 101.018           | Cortex Screw,                          | 2.7mm, SS, Length 18mm | 5 |       | 311.021           | Mini Chisels, Size: 10mm, SS                      |                                   |  | 1     |
| 101.020           | Cortex Screw,                          | 2.7mm, SS, Length 20mm | 5 |       | 311.022           | Mini Chisels, Size: 12mm, SS                      |                                   |  | 1     |
| 101.022           | Cortex Screw,                          | 2.7mm, SS, Length 22mm | 5 |       | 311.023           | Mini Osteotomes, Size: 4mm, SS                    |                                   |  | 1     |
| 101.024           | Cortex Screw,                          | 2.7mm, SS, Length 24mm | 5 |       | 311.024           | Mini Osteotomes, Size: 6mm, SS                    |                                   |  | 1     |
| 131.103           | Quarter Tubular Plate 2.7, SS 3 Holes  |                        | 2 |       | 311.025           | Mini Osteotomes, Size: 8mm, SS                    |                                   |  | 1     |
| 131.104           | Quarter Tubular Plate 2.7, SS, 4 Holes |                        | 2 |       | 311.026           | Mini Osteotomes, Size: 10mm, SS                   |                                   |  | 1     |
| 131.105           | Quarter Tubular Plate 2.7, SS, 5 Holes |                        | 2 |       | 311.027           | Mini Osteotomes, Size: 12mm, SS                   |                                   |  | 1     |
| 131.106           | Quarter Tubular Plate 2.7, SS, 6 Holes |                        | 1 |       | 311.028           | Mini Gouge, Size: 4mm, SS                         |                                   |  | 1     |
| 131.107           | Quarter Tubular Plate 2.7, SS, 7 Holes |                        | 1 |       | 311.029           | Mini Gouge, Size: 6mm, SS                         |                                   |  | 1     |
| 131.108           | Quarter Tubular Plate 2.7, SS, 8 Holes |                        | 1 |       | 311.030           | Mini Gouge, Size: 8mm, SS                         |                                   |  | 1     |
| 133.101           | L-Plate 2.7, SS, Left                  |                        | 2 |       | 311.031           | Mini Gouge, Size: 10mm, SS                        |                                   |  | 1     |
| 133.102           | L-Plate 2.7, SS, Right                 |                        | 2 |       | 311.032           | Mini Gouge, Size: 12mm, SS                        |                                   |  | 1     |
| 134.100           | T-Plate 2.7, SS                        |                        | 4 |       | 311.034           | Mini Hohmann Retractors, Size: 6mm, SS            |                                   |  | 1     |
| 121.103           | Straight Plate 2.0, SS, 3 Holes        |                        | 2 |       | 311.035           | Mini Hohmann Retractors, Size: 8mm, SS            |                                   |  | 1     |
| 121.104           | Straight Plate 2.0, SS, 4 Holes        |                        | 2 |       | 311.036           | Mini Hohmann Retractors, Size: 10mm, SS           |                                   |  | 1     |
| 121.105           | Straight Plate 2.0, SS, 5 Holes        |                        | 2 |       | 301.040           | Small Hexagonal Screw Driver, 2.5mm tip           |                                   |  | 1     |
| 121.106           | Straight Plate 2.0, SS, 6 Holes        |                        | 2 |       | 906.049           | Empty Container                                   |                                   |  |       |
| 122.101           | L-Plate 2.0, SS, Left                  |                        | 2 |       |                   |                                                   |                                   |  |       |
| 122.102           | L-Plate 2.0, SS, Right                 |                        | 2 |       |                   |                                                   |                                   |  |       |
| 123.100           | T-Plate 2.0, SS                        |                        | 4 |       |                   |                                                   |                                   |  |       |
| 74.03             | Condylar Plate 2.0, SS with Left Pin   |                        | 2 |       |                   |                                                   |                                   |  |       |
| 74.04             | Condylar Plate 2.0, SS with Right Pin  |                        | 2 |       |                   |                                                   |                                   |  |       |
| 74.01             | Condylar Plate 1.5, SS with Left Pin   |                        | 2 |       |                   |                                                   |                                   |  |       |
| 74.02             | Condylar Plate 1.5, SS with Right Pin  |                        | 2 |       |                   |                                                   |                                   |  |       |



**Mini Hexagonal Screw Driver (1.5mm Tip)  
with Screw Holding Sleeve**  
Code No.: 311.002

YouTube [www.TV.ortho.in/4016](http://www.TV.ortho.in/4016)



**Mini Screw Holding Forceps (for 1.5 & 2.0mm Screws)**  
Code No.: 311.011



**Mini Tap Sleeve**  
Code No.: 311.004

YouTube [www.TV.ortho.in/4566](http://www.TV.ortho.in/4566)



**Mini Periosteal Elevator**  
Code No.: 311.017

YouTube [www.TV.ortho.in/49](http://www.TV.ortho.in/49)



**Mini Depth Gouge**  
Code No.: 311.005

YouTube [www.TV.ortho.in/4567](http://www.TV.ortho.in/4567)



**Mini Drill Sleeve**  
Code No.: 311.003



**Mini Plate Bender (Pair)**  
Code No.: 311.010  
[YouTube www.TV.ortho.in/4149](https://www.youtube.com/watch?v=www.TV.ortho.in/4149)



**Mini Osteotomes**  
Sizes: 4, 6, 8, 10 & 12mm  
Code No.: 311.023 to 027  
[YouTube www.TV.ortho.in/5894](https://www.youtube.com/watch?v=www.TV.ortho.in/5894)



**Mini Gouge**  
Sizes: 4, 6, 8, 10 & 12mm  
Code No.: 311.028 to 032  
[YouTube www.TV.ortho.in/16492](https://www.youtube.com/watch?v=www.TV.ortho.in/16492)



**Mini Self Retaining Retractor for Fingers**  
Code No.: 311.006  
[YouTube www.TV.ortho.in/4568](https://www.youtube.com/watch?v=www.TV.ortho.in/4568)



**Mini Self Retaining Retractor for Metacarpals**  
Code No.: 311.007  
[YouTube www.TV.ortho.in/4146](https://www.youtube.com/watch?v=www.TV.ortho.in/4146)



**Mini Plate Holding Forceps**  
Code No.: 311.013  
[YouTube www.TV.ortho.in/4152](https://www.youtube.com/watch?v=www.TV.ortho.in/4152)



**Mini Reduction Forceps (Serrated)**  
Code No.: 311.015  
[YouTube www.TV.ortho.in/4154](https://www.youtube.com/watch?v=www.TV.ortho.in/4154)



**Mini Neutral & Loaded Drill Guide**  
Code No.: 311.009  
[YouTube www.TV.ortho.in/4148](https://www.youtube.com/watch?v=www.TV.ortho.in/4148)



**Mini Chisels**  
Sizes: 4, 6, 8, 10 & 12mm  
Code No.: 311.018 to 022  
[YouTube www.TV.ortho.in/5890](https://www.youtube.com/watch?v=www.TV.ortho.in/5890)



**Mini Hohmann Retractors**  
Sizes: 6, 8, 10mm  
Code No.: 311.034 to 036  
[YouTube www.TV.ortho.in/16494](https://www.youtube.com/watch?v=www.TV.ortho.in/16494)



**Hollow Mill**  
for 1.5mm & 2.0mm Screws  
Code No.: 311.012



**Mini Self Retaining Retractor for Metacarpals (with Hinge)**  
Code No.: 311.008  
[YouTube www.TV.ortho.in/4147](https://www.youtube.com/watch?v=www.TV.ortho.in/4147)



**Mini Plate & Bone Holding Forceps**  
Code No.: 311.014  
[YouTube www.TV.ortho.in/4153](https://www.youtube.com/watch?v=www.TV.ortho.in/4153)



**Mini Reduction Forceps (Pointed)**  
Code No.: 311.016  
[YouTube www.TV.ortho.in/4155](https://www.youtube.com/watch?v=www.TV.ortho.in/4155)