



□ "S" SERIES - SAFETY LOCKING IMPLANTS




S-SERIES LOCKING IMPLANTS

...SAFETY LOCK PLATES



FILOS-PROXIMAL HUMERUS
SAFETY LOCK PLATE 3.5 - STANDARD

No. of Holes	Length in mm	Code TITANIUM
3	90	sTi-SF-642.003
5	114	sTi-SF-642.005



FILOS-PROXIMAL HUMERUS
SAFETY LOCK PLATE 3.5 - LONG

No. of Holes	Length in mm	Code TITANIUM
5	142	sTi-SF-642.205
6	160	sTi-SF-642.206
8	196	sTi-SF-642.208
10	232	sTi-SF-642.210
12	268	sTi-SF-642.212



SAFETY LOCK 'T' PLATE 3.5, OBLIQUE ANGLED

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
3	50	sTi-SF-142.103L	sTi-SF-142.103R
4	57	sTi-SF-142.104L	sTi-SF-142.104R
5	67	sTi-SF-142.105L	sTi-SF-142.105R
6	77	sTi-SF-142.106L	sTi-SF-142.106R
7	87	sTi-SF-142.107L	sTi-SF-142.107R
8	97	sTi-SF-142.108L	sTi-SF-142.108R



PROXIMAL HUMERUS
SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	Code TITANIUM
5	85	sTi-SF-640.005
8	121	sTi-SF-640.008



SAFETY LOCK 'T' PLATE 3.5 - RIGHT ANGLED

No. of Holes	Code TITANIUM
3	sTi-SF-141.103
4	sTi-SF-141.104
5	sTi-SF-141.105
6	sTi-SF-141.106
7	sTi-SF-141.107
8	sTi-SF-141.108



LC-DCP SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	Code TITANIUM
4	59 mm	sTi-SF-147.104
5	72 mm	sTi-SF-147.105
6	85 mm	sTi-SF-147.106
7	98 mm	sTi-SF-147.107
8	111 mm	sTi-SF-147.108
9	124 mm	sTi-SF-147.109
10	137 mm	sTi-SF-147.110
11	150 mm	sTi-SF-147.111
12	163 mm	sTi-SF-147.112



BLOCK FOR FILOS-PROXIMAL HUMERUS
SAFETY LOCK PLATE
Code No.: BL-642




BLOCK FOR SAFETY LOCKING
DISTAL TIBIAL PLATE
Code No.: BL-646



MEDIAL DISTAL TIBIA SAFETY LOCK PLATE 3.5 (WITH TAB)

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
4	116	sTi-SF-645.04L	sTi-SF-645.04R
6	142	sTi-SF-645.06L	sTi-SF-645.06R
8	168	sTi-SF-645.08L	sTi-SF-645.08R
10	194	sTi-SF-645.10L	sTi-SF-645.10R
12	220	sTi-SF-645.12L	sTi-SF-645.12R
14	246	sTi-SF-645.14L	sTi-SF-645.14R



MEDIAL DISTAL TIBIAL SAFETY LOCK PLATE 3.5 WITHOUT TAB

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
4	116	sTi-SF-646.04L	sTi-SF-646.04R
6	142	sTi-SF-646.06L	sTi-SF-646.06R
8	168	sTi-SF-646.08L	sTi-SF-646.08R
10	194	sTi-SF-646.10L	sTi-SF-646.10R
12	220	sTi-SF-646.12L	sTi-SF-646.12R
14	246	sTi-SF-646.14L	sTi-SF-646.14R



METAPHYSEAL SAFETY LOCK PLATE 3.5/4.5/5.0 - BROAD

No. of Holes	Length in mm	Code TITANIUM
5+3	118	sTi-SF-615.03
5+4	136	sTi-SF-615.04
5+5	154	sTi-SF-615.05
5+6	172	sTi-SF-615.06
5+7	190	sTi-SF-615.07
5+8	208	sTi-SF-615.08
5+9	226	sTi-SF-615.09
5+11	262	sTi-SF-615.11
5+13	298	sTi-SF-615.13
5+15	334	sTi-SF-615.15



CLAVICLE HOOK SAFETY LOCK PLATE 3.5

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
4	12	sTi-SF-609.01L	sTi-SF-609.01R
4	15	sTi-SF-609.03L	sTi-SF-609.03R
4	18	sTi-SF-609.05L	sTi-SF-609.05R
5	12	sTi-SF-609.07L	sTi-SF-609.07R
5	15	sTi-SF-609.09L	sTi-SF-609.09R
5	18	sTi-SF-609.11L	sTi-SF-609.11R
6	15	sTi-SF-609.13L	sTi-SF-609.13R
6	18	sTi-SF-609.15L	sTi-SF-609.15R
7	15	sTi-SF-609.17L	sTi-SF-609.17R
7	18	sTi-SF-609.19L	sTi-SF-609.19R



LATERAL DISTAL FIBULA SAFETY LOCK PLATE 2.7/3.5

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
3	73	sTi-SF-652.03L	sTi-SF-652.03R
4	86	sTi-SF-652.04L	sTi-SF-652.04R
5	99	sTi-SF-652.05L	sTi-SF-652.05R
6	112	sTi-SF-652.06L	sTi-SF-652.06R
7	125	sTi-SF-652.07L	sTi-SF-652.07R
8	138	sTi-SF-652.08L	sTi-SF-652.08R



RECONSTRUCTION SAFETY LOCK PLATE 3.5, ROUND HOLES - STRAIGHT

No. of Holes	Length in mm	Code TITANIUM
5	65	sTi-SF-444.105
6	78	sTi-SF-444.106
7	91	sTi-SF-444.107
8	104	sTi-SF-444.108
9	117	sTi-SF-444.109
10	130	sTi-SF-444.110
12	156	sTi-SF-444.112
14	182	sTi-SF-444.114
16	208	sTi-SF-444.116
18	234	sTi-SF-444.118
20	260	sTi-SF-444.120
22	286	sTi-SF-444.122



RECONSTRUCTION SAFETY LOCK PLATE 3.5 STRAIGHT

No. of Holes	Length in mm	Code TITANIUM
5	70	sTi-SF-144.105
6	84	sTi-SF-144.106
7	98	sTi-SF-144.107
8	112	sTi-SF-144.108
9	126	sTi-SF-144.109
10	140	sTi-SF-144.110
12	168	sTi-SF-144.112
14	196	sTi-SF-144.114



**SAFETY LOCK SCREW 3.5MM
(SELF TAPPING)**

Length in mm	Code TITANIUM
10	sTi-SF-102.210
12	sTi-SF-102.212
14	sTi-SF-102.214
16	sTi-SF-102.216
18	sTi-SF-102.218
20	sTi-SF-102.220
22	sTi-SF-102.222
24	sTi-SF-102.224
26	sTi-SF-102.226
28	sTi-SF-102.228
30	sTi-SF-102.230
35	sTi-SF-102.235
40	sTi-SF-102.240
45	sTi-SF-102.245
50	sTi-SF-102.250
55	sTi-SF-102.255
60	sTi-SF-102.260
65	sTi-SF-102.265
70	sTi-SF-102.270
75	sTi-SF-102.275
80	sTi-SF-102.280
85	sTi-SF-102.285
90	sTi-SF-102.290



**SAFETY LOCK SCREW 2.7MM
(SELF TAPPING)**

Length in mm	Code TITANIUM
20	sTi-SF-101.420
22	sTi-SF-101.422
24	sTi-SF-101.424
26	sTi-SF-101.426
28	sTi-SF-101.428
30	sTi-SF-101.430
32	sTi-SF-101.432
34	sTi-SF-101.434
36	sTi-SF-101.436
38	sTi-SF-101.438
40	sTi-SF-101.440
42	sTi-SF-101.442
44	sTi-SF-101.444
46	sTi-SF-101.446
48	sTi-SF-101.448
50	sTi-SF-101.450
52	sTi-SF-101.452
54	sTi-SF-101.454
56	sTi-SF-101.456
58	sTi-SF-101.458
60	sTi-SF-101.460



**CORTEX SCREW 3.5MM
(SELF TAPPING)**

Length in mm	Code TITANIUM
10	sTi-102.210
12	sTi-102.212
14	sTi-102.214
16	sTi-102.216
18	sTi-102.218
20	sTi-102.220
22	sTi-102.222
24	sTi-102.224
26	sTi-102.226
28	sTi-102.228
30	sTi-102.230
32	sTi-102.232
34	sTi-102.234
36	sTi-102.236
38	sTi-102.238
40	sTi-102.240
45	sTi-102.245
50	sTi-102.250
55	sTi-102.255
60	sTi-102.260
65	sTi-102.265
70	sTi-102.270



SPACER

Size in mm	Code TITANIUM
3.5	sTi-115.35
5.0	sTi-115.50



C.C. WASHER

Size in mm	Code TITANIUM
4.0	sTi-115.40
6.5	sTi-115.65



**WIRE MOUNT - 5.0MM
TITANIUM**

Code No.: sTi-115.45



**CANCELLOUS BONE SCREW 4.0MM
SHORT THREAD**

Length in mm	Code TITANIUM
10	sTi-104.010
12	sTi-104.012
14	sTi-104.014
16	sTi-104.016
18	sTi-104.018
20	sTi-104.020
22	sTi-104.022
24	sTi-104.024
26	sTi-104.026
28	sTi-104.028
30	sTi-104.030
32	sTi-104.032
35	sTi-104.035
40	sTi-104.040
45	sTi-104.045
50	sTi-104.050
55	sTi-104.055
60	sTi-104.060
65	sTi-104.065
70	sTi-104.070
75	sTi-104.075



**CANCELLOUS BONE SCREW 4.0MM
FULL THREAD**

Length in mm	Code TITANIUM
10	sTi-105.010
12	sTi-105.012
14	sTi-105.014
16	sTi-105.016
18	sTi-105.018
20	sTi-105.020
22	sTi-105.022
24	sTi-105.024
26	sTi-105.026
28	sTi-105.028
30	sTi-105.030
32	sTi-105.032
35	sTi-105.035
40	sTi-105.040
45	sTi-105.045
50	sTi-105.050
55	sTi-105.055
60	sTi-105.060
65	sTi-105.065
70	sTi-105.070
75	sTi-105.075



CANCELLOUS CANNULATED SCREW 4.0MM - SHORT THREAD

Length in mm	Code TITANIUM
10	sTi-116.010
12	sTi-116.012
14	sTi-116.014
16	sTi-116.016
18	sTi-116.018
20	sTi-116.020
22	sTi-116.022
24	sTi-116.024
26	sTi-116.026
28	sTi-116.028
30	sTi-116.030
32	sTi-116.032
34	sTi-116.034
36	sTi-116.036
38	sTi-116.038
40	sTi-116.040
42	sTi-116.042
44	sTi-116.044
46	sTi-116.046
48	sTi-116.048
50	sTi-116.050
55	sTi-116.055
60	sTi-116.060

□ SMALL FRAGMENT SAFETY LOCKING INSTRUMENTS SET

sAS432 SMALL FRAGMENT SAFETY LOCKING INSTRUMENTS SET

	<i>Set Consisting of:</i>	<i>Units</i>
sAS432.003	Drill Bit, 2.5 mm	1
sAS432.001	Drill Bit, 2.8 mm	2
sAS432.011	Tap, 3.5 mm	1
sAS432.015	Screw Driver Shaft, 3.5 mm	1
sAS432.005	Drill Sleeve, 2.8 mm	2
sAS432.013	T Handle With Quick Coupling	1
sAS432.007	Double Drill Sleeve 2.5/3.5 mm	1
sAS432.023	Depth Gauge	1
sAS432.019	Screw Driver, 2.5 mm Tip	1
sAS432.009	Torque Limiting Screw Driver, 2.5 mm Tip	1
sAS432.029	Empty Container	



Drill Bit, 2.5mm
Code No.: sAS432.003



Drill Bit, 2.8mm
Code No.: sAS432.001



Tap, 3.5mm
Code No.: sAS432.011



T Handle With Quick Coupling
Code No.: sAS432.013



Drill Sleeve, 2.8mm
Code No.: sAS432.005



Screw Driver Shaft, 3.5mm
Code No.: sAS432.015



Double Drill Sleeve 2.5/3.5mm
Code No.: sAS432.007



Depth Gauge
Code No.: sAS432.023



Screw Driver, 2.5mm Tip
Code No.: sAS432.019



Torque Limiting Screw Driver, 2.5mm Tip
Code No.: sAS432.009



NARROW LC-DCP SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Length in mm	Code TITANIUM
4	80 mm	sTi-SF-150.104
5	98 mm	sTi-SF-150.105
6	116 mm	sTi-SF-150.106
7	134 mm	sTi-SF-150.107
8	152 mm	sTi-SF-150.108
9	170 mm	sTi-SF-150.109
10	188 mm	sTi-SF-150.110
11	206 mm	sTi-SF-150.111
12	224 mm	sTi-SF-150.112
14	260 mm	sTi-SF-150.114
16	296 mm	sTi-SF-150.116



BROAD LC-DCP SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Length in mm	Code TITANIUM
5	98 mm	sTi-SF-151.105
6	116 mm	sTi-SF-151.106
7	134 mm	sTi-SF-151.107
8	152 mm	sTi-SF-151.108
9	170 mm	sTi-SF-151.109
10	188 mm	sTi-SF-151.110
11	206 mm	sTi-SF-151.111
12	224 mm	sTi-SF-151.112
13	242 mm	sTi-SF-151.113
14	260 mm	sTi-SF-151.114
15	278 mm	sTi-SF-151.115
16	296 mm	sTi-SF-151.116
18	332 mm	sTi-SF-151.117



MEDIAL PROXIMAL TIBIAL SAFETY LOCK PLATE 4.5

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
4	106	sTi-SF-653.004L	sTi-SF-653.004R
6	142	sTi-SF-653.006L	sTi-SF-653.006R
8	178	sTi-SF-653.008L	sTi-SF-653.008R
10	214	sTi-SF-653.010L	sTi-SF-653.010R
12	250	sTi-SF-653.012L	sTi-SF-653.012R
14	286	sTi-SF-653.014L	sTi-SF-653.014R



PROXIMAL FEMORAL SAFETY LOCK PLATE 4.5/5.0/7.3

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
4	175	sTi-SF-656.004L	sTi-SF-656.004R
6	211	sTi-SF-656.006L	sTi-SF-656.006R
8	247	sTi-SF-656.008L	sTi-SF-656.008R
10	283	sTi-SF-656.010L	sTi-SF-656.010R
12	319	sTi-SF-656.012L	sTi-SF-656.012R
14	355	sTi-SF-656.014L	sTi-SF-656.014R
16	391	sTi-SF-656.016L	sTi-SF-656.016R




SAFETY LOCK 'T' PLATE 4.5/5.0

No. of Holes	Length in mm	Code TITANIUM
3	67	sTi-SF-159.103
4	83	sTi-SF-159.104
5	99	sTi-SF-159.105
6	115	sTi-SF-159.106
7	131	sTi-SF-159.107
8	147	sTi-SF-159.108
9	163	sTi-SF-159.109
10	179	sTi-SF-159.110
12	195	sTi-SF-159.112




T BUTTRESS SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Code TITANIUM
4	sTi-SF-160.104
5	sTi-SF-160.105
6	sTi-SF-160.106
7	sTi-SF-160.107
8	sTi-SF-160.108
9	sTi-SF-160.109
10	sTi-SF-160.110
12	sTi-SF-160.112



PROXIMAL LATERAL TIBIAL SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
5	140	sTi-SF-165.105L	sTi-SF-165.005R
7	180	sTi-SF-165.107L	sTi-SF-165.007R
9	220	sTi-SF-165.109L	sTi-SF-165.009R
11	260	sTi-SF-165.111L	sTi-SF-165.011R
13	300	sTi-SF-165.113L	sTi-SF-165.013R



DISTAL FEMORAL SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Length mm	Code - TITANIUM	
		Left	Right
5	156	sTi-SF-635.005L	sTi-SF-635.005R
7	196	sTi-SF-635.007L	sTi-SF-635.007R
9	236	sTi-SF-635.009L	sTi-SF-635.009R
11	276	sTi-SF-635.011L	sTi-SF-635.011R
13	316	sTi-SF-635.013L	sTi-SF-635.013R



BLOCK FOR PROXIMAL LATERAL TIBIAL SAFETY LOCK PLATE

Code No.: BL-165



BLOCK FOR DISTAL FEMORAL SAFETY LOCK PLATE

Code No.: BL-635



L. BUTTRESS SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Code - TITANIUM	
	Left	Right
3	sTi-SF-161.133L	sTi-SF-161.103R
4	sTi-SF-161.134L	sTi-SF-161.104R
5	sTi-SF-161.135L	sTi-SF-161.105R
6	sTi-SF-161.136L	sTi-SF-161.106R
7	sTi-SF-161.137L	sTi-SF-161.107R
8	sTi-SF-161.138L	sTi-SF-161.108R
9	sTi-SF-161.139L	sTi-SF-161.109R
10	sTi-SF-161.140L	sTi-SF-161.110R
12	sTi-SF-161.142L	sTi-SF-161.112R



SAFETY LOCK SCREW 5.0MM
(SELF TAPPING)

Length in mm	Code TITANIUM
14	sTi-SF-500.014
16	sTi-SF-500.016
18	sTi-SF-500.018
20	sTi-SF-500.020
22	sTi-SF-500.022
24	sTi-SF-500.024
26	sTi-SF-500.026
28	sTi-SF-500.028
30	sTi-SF-500.030
32	sTi-SF-500.032
34	sTi-SF-500.034
36	sTi-SF-500.036
38	sTi-SF-500.038
40	sTi-SF-500.040
42	sTi-SF-500.042
44	sTi-SF-500.044
46	sTi-SF-500.046
48	sTi-SF-500.048
50	sTi-SF-500.050
55	sTi-SF-500.055
60	sTi-SF-500.060
65	sTi-SF-500.065
70	sTi-SF-500.070
75	sTi-SF-500.075
80	sTi-SF-500.080
85	sTi-SF-500.085
90	sTi-SF-500.090
95	sTi-SF-500.095
100	sTi-SF-500.100
105	sTi-SF-500.105
110	sTi-SF-500.110



**CANNULATED CANCELLOUS
SAFETY LOCK SCREW 7.3MM**

Length in mm	Code TITANIUM
40	sTi-SF-505.740
45	sTi-SF-505.745
50	sTi-SF-505.750
55	sTi-SF-505.755
60	sTi-SF-505.760
65	sTi-SF-505.765
70	sTi-SF-505.770
75	sTi-SF-505.775
80	sTi-SF-505.780
85	sTi-SF-505.785
90	sTi-SF-505.790
95	sTi-SF-505.795
100	sTi-SF-505.800
105	sTi-SF-505.805
110	sTi-SF-505.810
115	sTi-SF-505.815



CORTEX SCREW 4.5MM
(SELF TAPPING)

Length in mm	Code TITANIUM
14	sTi-106.214
16	sTi-106.216
18	sTi-106.218
20	sTi-106.220
22	sTi-106.222
24	sTi-106.224
26	sTi-106.226
28	sTi-106.228
30	sTi-106.230
32	sTi-106.232
34	sTi-106.234
36	sTi-106.236
38	sTi-106.238
40	sTi-106.240
42	sTi-106.242
44	sTi-106.244
46	sTi-106.246
48	sTi-106.248
50	sTi-106.250
52	sTi-106.252
54	sTi-106.254
56	sTi-106.256
58	sTi-106.258
60	sTi-106.260
62	sTi-106.262
64	sTi-106.264
66	sTi-106.266
68	sTi-106.268
70	sTi-106.270
75	sTi-106.275
80	sTi-106.280
85	sTi-106.285
90	sTi-106.290



**CANCELLOUS BONE SCREW 6.5MM
FULL THREAD**

Length in mm	Code TITANIUM
25	sTi-110.025
30	sTi-110.030
35	sTi-110.035
40	sTi-110.040
45	sTi-110.045
50	sTi-110.050
55	sTi-110.055
60	sTi-110.060
65	sTi-110.065
70	sTi-110.070
75	sTi-110.075
80	sTi-110.080
85	sTi-110.085
90	sTi-110.090
95	sTi-110.095
100	sTi-110.100
105	sTi-110.105
110	sTi-110.110



**CANCELLOUS BONE SCREW 6.5MM
16MM THREAD**

Length in mm	Code TITANIUM
25	sTi-108.025
30	sTi-108.030
35	sTi-108.035
40	sTi-108.040
45	sTi-108.045
50	sTi-108.050
55	sTi-108.055
60	sTi-108.060
65	sTi-108.065
70	sTi-108.070
75	sTi-108.075
80	sTi-108.080
85	sTi-108.085
90	sTi-108.090
95	sTi-108.095
100	sTi-108.100
105	sTi-108.105
110	sTi-108.110



**CANCELLOUS BONE SCREW 6.5MM
32MM THREAD**

Length in mm	Code TITANIUM
40	sTi-109.040
45	sTi-109.045
50	sTi-109.050
55	sTi-109.055
60	sTi-109.060
65	sTi-109.065
70	sTi-109.070
75	sTi-109.075
80	sTi-109.080
85	sTi-109.085
90	sTi-109.090
95	sTi-109.095
100	sTi-109.100
105	sTi-109.105
110	sTi-109.110



CANCELLOUS CANNULATED SCREW
6.5MM - 16MM THREAD

Length in mm	Code TITANIUM
30	sTi-455.030
35	sTi-455.035
40	sTi-455.040
45	sTi-455.045
50	sTi-455.050
55	sTi-455.055
60	sTi-455.060
65	sTi-455.065
70	sTi-455.070
75	sTi-455.075
80	sTi-455.080
85	sTi-455.085
90	sTi-455.090
95	sTi-455.095
100	sTi-455.100
105	sTi-455.105
110	sTi-455.110



CANCELLOUS CANNULATED SCREW
6.5MM - 32MM THREAD

Length in mm	Code TITANIUM
40	sTi-455.240
45	sTi-455.245
50	sTi-455.250
55	sTi-455.255
60	sTi-455.260
65	sTi-455.265
70	sTi-455.270
75	sTi-455.275
80	sTi-455.280
85	sTi-455.285
90	sTi-455.290
95	sTi-455.295
100	sTi-455.300
105	sTi-455.305
110	sTi-455.310



CANCELLOUS CANNULATED SCREW
7.0MM- 16MM THREAD

Length in mm	Code TITANIUM
50	sTi-111.050
55	sTi-111.055
60	sTi-111.060
65	sTi-111.065
70	sTi-111.070
75	sTi-111.075
80	sTi-111.080
85	sTi-111.085
90	sTi-111.090
95	sTi-111.095
100	sTi-111.100
105	sTi-111.105
110	sTi-111.110



CANCELLOUS CANNULATED SCREW
7.0MM - 32MM THREAD

Length in mm	Code TITANIUM
50	sTi-112.050
55	sTi-112.055
60	sTi-112.060
65	sTi-112.065
70	sTi-112.070
75	sTi-112.075
80	sTi-112.080
85	sTi-112.085
90	sTi-112.090
95	sTi-112.095
100	sTi-112.100
105	sTi-112.105
110	sTi-112.110



**DHS SAFETY LOCK PLATE 135°
38MM BARREL**

No. of Holes	Code TITANIUM
4	sTi-SF-166.024
5	sTi-SF-166.025
6	sTi-SF-166.026
7	sTi-SF-166.027
8	sTi-SF-166.028
9	sTi-SF-166.029
10	sTi-SF-166.030
12	sTi-SF-166.032
14	sTi-SF-166.034
16	sTi-SF-166.036



**DHS SAFETY LOCK PLATE 135°
SHORT BARREL**

No. of Holes	Code TITANIUM
4	sTi-SF-166.084
5	sTi-SF-166.085
6	sTi-SF-166.086
7	sTi-SF-166.087
8	sTi-SF-166.088
9	sTi-SF-166.089
10	sTi-SF-166.090
12	sTi-SF-166.092
14	sTi-SF-166.094
16	sTi-SF-166.096



DCS SAFETY LOCK PLATE 95°

No. of Holes	Code TITANIUM
4	sTi-SF-167.004
5	sTi-SF-167.005
6	sTi-SF-167.006
7	sTi-SF-167.007
8	sTi-SF-167.008
9	sTi-SF-167.009
10	sTi-SF-167.010
12	sTi-SF-167.012
14	sTi-SF-167.014
16	sTi-SF-167.016



**STANDARD LAG SCREW 12.5MM
(WITH COMPRESSION SCREW)**

Length in mm	Code TITANIUM
50	sTi-168.050
55	sTi-168.055
60	sTi-168.060
65	sTi-168.065
70	sTi-168.070
75	sTi-168.075
80	sTi-168.080
85	sTi-168.085
90	sTi-168.090
95	sTi-168.095
100	sTi-168.100
105	sTi-168.105
110	sTi-168.110
115	sTi-168.115

□ LARGE FRAGMENT SAFETY LOCKING INSTRUMENTS SET

sAS435 LARGE FRAGMENT SAFETY LOCKING INSTRUMENTS SET

Set Consisting of:		Units
sAS435.004	Drill Bit, 3.2 mm	1
sAS435.002	Drill Bit, 4.3 mm	2
sAS435.012	Tap, 4.5 mm	1
sAS435.016	Screw Driver Shaft, 4.5 mm	1
sAS435.006	Drill Sleeve, 4.3 mm	2
sAS435.014	T Handle With Quick Coupling	1
sAS435.008	Double Drill Sleeve 3.2/4.5 mm	1
sAS435.024	Depth Gauge	1
sAS435.020	Screw Driver, 3.5 mm Tip	1
sAS435.010	Torque Limiting Screw Driver, 3.5 mm Tip	1
sAS435.030	Empty Container	



Drill Bit, 3.2mm
Code No.: sAS435.004



Drill Bit, 4.3mm
Code No.: sAS435.002



Tap, 4.5mm
Code No.: sAS435.012



T Handle With Quick Coupling
Code No.: sAS435.014



Drill Sleeve, 4.3mm
Code No.: sAS435.006



Screw Driver Shaft, 4.5mm
Code No.: sAS435.016



Double Drill Sleeve 3.2/4.5mm
Code No.: sAS435.008



Depth Gauge
Code No.: sAS435.024



Screw Driver, 3.5mm Tip
Code No.: sAS435.020



Torque Limiting Screw Driver, 3.5mm Tip
Code No.: sAS435.010